

Gas Scrubber 가스 스크러버

Cat. No.	Model	Description	EA Price ₩
G01-132-050	1830-10	3 L Gas Scrubber Set	11,267,700
G01-132-055	1830-50	20 L Gas Scrubber Set	22,462,900


- 여러 종류의 Gas를 정화하는데 사용됨
- 3 L와 20 L, 2종류의 Reservoir 사이즈
- 3 L Set는 Hood내에서 사용하기에 적합한 크기
- Gas Flow (1830-10) Max. 3 L/min
(1830-50) Max. 8 L/min
- 구성 : Reservoir, Absorption Column w/PTFE Disc, Glass Rasching Ring Packing Exhaust Gas Outlet Column, Top Shower Head Adapter, Liquid Transfer Pump, 3-P Clamp, 3ea Gasket, Condenser, Gas Inlet Adapter, PTFE Tubing, 5ea PFA Fitting
- Stand는 별도
- Chemglass

Optional Stand

Cat. No.	Model	Description	적용 Model	EA Price ₩
G01-132-091	1949-100	Support Stand, Bench Top	1830-10	5,858,700
G01-132-095	1951-100	Support Stand, Work In Hood	1830-50	9,086,700
G01-132-099	1830-60			9,259,500



1830-50



1949-100

Gas Manifold 가스 메니폴드

Cat. No.	Model	Cock	전체 길이	EA Price ₩
G01-132-145	1814-01	Glass, Bore 2 mm	365 mm	684,700
G01-132-151	1816-01	Teflon, Bore 2 mm		666,000

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Capillary Tubing Size : od8×id 2 mm
- Chemglass



Glass Cock



PTFE Cock

Midget Impinger 임핀저

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 소용량의 집진 장치, EPA Method 6
- Bottle과 Stopper Joint Size : 24/40
- Chemglass

With 1 mm Impinger tip Orifice

Cat. No.	Model	용량	Tubulation	EA Price ₩
G01-132-157	1820-01 ①	30 ml	Plain	147,100
G01-132-163	1820-02 ②		12/5	298,300


- 0.09 ~ 0.11 CFM at 12" H₂O Vacuum

With Fritted Gas Dispersion Tube

Cat. No.	Model	용량	Tubulation	EA Price ₩
G01-132-169	1821-01 ③	30 ml	Plain	187,500



1820-01



1820-02



1821-01