

## Gas Dispersion Tubes

### 가스 분산튜브

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Porous Fritted Glass
- Chemglass

#### Tube od6 mm, Overall Height 110 mm ①

Cat. No.	Model	Porosity	Fritter Size	EA Price ₩
G01-132-211	207-01	F, 4 ~ 5.5 um	Tube Type, od7×h12 mm	32,000
G01-132-217	207-02	M, 10 ~ 16 um		32,000
G01-132-223	207-03	C, 40 ~ 60 um		32,000
G01-132-229	207-04	XC, 170 ~ 220 um		34,500

#### Tube od8 mm, Overall Height 250 mm ②

Cat. No.	Model	Porosity	Fritter Size	EA Price ₩
G01-132-235	203-01	F, 4 ~ 5.5 um	Cylinder Type, od12×h20 mm	45,500
G01-132-241	203-02	M, 10 ~ 16 um		45,500
G01-132-247	203-03	C, 40 ~ 60 um		45,500
G01-132-253	203-04	XC, 170 ~ 220 um		45,500

#### Tube od8 mm, Overall Height 250 mm ③

Cat. No.	Model	Porosity (um)	Filter od	전체 Size (mm)	EA Price ₩
G01-132-265	204-06	M, 10 ~ 16	20 mm	od8×h115	59,600
G01-132-261	204-07	C, 40 ~ 60			59,600
G01-132-279	204-10	M, 10 ~ 16	30 mm	od8×h190	61,800
G01-132-285	204-11	C, 40 ~ 60			61,800



## Bubbler, Gas Distribution Tube

### 버블러, 가스 분산 튜브, LukeGL®

- LukeGL®

#### Cylindrical Type, Tube Ø8×L210 mm ①

Cat. No.	Model	Filter Plug	Pore Size	EA Price ₩
B11-10-240	7208020	Ø9×20 mm	P0, 160 ~ 250 um	65,800
B11-10-250	7208120		P1, 100 ~ 160 um	58,800
B11-10-260	7208220		P2, 40 ~ 100 um	58,800
B11-10-270	7208320		P3, 13 ~ 40 um	65,800
B11-10-280	7208420		P4, 10 ~ 16 um	75,400
B11-10-290	7209025	Ø12×25 mm	P0, 160 ~ 250 um	70,600
B11-10-300	7209125		P1, 100 ~ 160 um	65,800
B11-10-310	7209225		P2, 40 ~ 100 um	70,600
B11-10-320	7209325		P3, 13 ~ 40 um	70,600
B11-10-330	7209425		P4, 10 ~ 16 um	76,700

#### Filter Plug Type, Ø9×L210 mm

Cat. No.	Model	Filter Plug	Pore Size	EA Price ₩
B11-10-340	7211025	Ø22 mm 상향식 ②	P0, 160 ~ 250 um	20,200
B11-10-350	7211125		P1, 100 ~ 160 um	20,200
B11-10-360	7211225		P2, 40 ~ 100 um	72,400
B11-10-370	7212025	Ø25 mm 하향식 ③	P0, 160 ~ 250 um	20,600
B11-10-380	7212125		P1, 100 ~ 160 um	24,500
B11-10-390	7212225		P2, 40 ~ 100 um	72,400

