

Gas Drying Unit, Indicating 가스 건조관

Cat. No.	Model	Description (Type)	EA Price ₩
D08-55-411	26800	Hose Barbs for Flexible Tubing ①	206,700
D08-55-417	26840	1/8" SUS Female NPT ②	346,000
D08-55-423	26809	Mounting Clip (Accessory) ④	16,000
D08-55-471	26805	Filter Disc for Gas Drying Unit ⑤	2,200

- 재질 : Polycarbonate, Stainless Steel Support 와 Spring
- Size : Ø67×h290 mm
- Max. Working Pressure : 90 psig, 200 L/hr
- Desiccant : 8 mesh Indicating DRIERITE (23001, 23005)
- Drierite



Gas Purifier, Indicating 가스 정제관, 가스정제 및 건조겸용

Cat. No.	Model	Description (Type)	EA Price ₩
D08-55-455	27068	Gas Purifier, up to 300 L/hr ③	443,400
D08-55-423	26809	Mounting Clip (Accessory) ④	16,000
D08-55-471	26805	Filter Disc for Gas Drying Unit ⑤	2,200
D08-55-467	27069	Refill Kit	105,100

- Size : Ø66.7×h289 mm
- Max. Working Pressure : 100 psig, 300 L/h
- PC케이스 및 Cover와 316SS 40 Micro Outlet Frit 포함
- Connection : 1/8" Stainless Steel Male Tube Fitting
- Desiccant : 250 cc 8 mesh Indicating DRIEREIT & 250 cc 5A Molecular Sieves
- 가스의 정제 (5A Molecular까지) 및 습기 제거용 (25 g까지)
- Drierite



Gas Drying Tube, Indicating 가스 건조 튜브

Cat. No.	Model	Description	EA Price ₩
D08-55-522	26930	Gas Drying Tube	28,800

- Size : Ø19×L200 mm (Tube Length : 150 mm)
- Connection : Hose Barbs for Flexible Tubing
- Desiccant : 30 g / 8 mesh Indicating DRIERITE (23001, 23005)
- 재질
(Tube) : PC (Polycarbonate)
(Cap) : Molded PP (Polypropylene)
- Drierite



Drierite Absorbent 흡습 건조제, Indicating

Cat. No.	Model	용량	Mesh	EA Price ₩	CS Price ₩
D08-55-354	23001	450 g	8	50,800	12ea 569,800
D08-55-360	23005	2.25 kg		248,800	4ea 929,900
D08-55-371	24001	450 g	10 ~ 20	51,700	12ea 579,900
D08-55-373	24005	2.25 kg		253,000	4ea 945,800

- 수분 흡수율 : 10 ~ 14% in Their Weight
- 210 ~ 230°C 에서 1시간 가량 건조 시키면 재사용이 가능함
- 건조제 : 수분 흡수 물질 (색상 변화 : Blue → Pink)
- Contents : < 98% CaSO₄, >2% CoCl₂
- Drierite

