


## SPure™ Inline Gas Filter, Plastic Body 플라스틱 케이스 가스 정제관

- General Purpose용 Moisture, Oil, Dust & Hydrocarbon 트랩
- 재질 : PC or Acrylic(400cc only) Body, Aluminum Endcap
- Compressed Fitting : Brass or SS(Optional), Push-to-connect Fitting
- 고압의 Compressed Air Line이나 GC Carrier Gas와 같이 사용하기에 부적합
- Hazardous, Flammable, or Reactive Gas와 같이 사용하기에 부적합
- Cobalt-free Moisture Indicator Signal : Green to Blue, at 5% Relative Humidity
- Chromres



Capacity	88cc	116cc	240cc	400cc
Housing Material	Polycarbonate			Acrylic
Size, Compression Fitting	∅38 x h263 mm	∅38 x h318 mm	∅56 x h285 mm	∅56 x h419 mm
Size, Push-to-connect Fitting	∅38 x h243 mm	∅38 x h298 mm	∅56 x h265 mm	∅56 x h399 mm
Max. Pressure	6.89 bar / 100 psi			

### SPure™ H<sub>2</sub>O Purifier, 플라스틱 케이스 수분 트랩 Indicator : Green to Blue (at 5% Relative Humidity)

Cat. No.	Model	용량	Trap 용량	Fitting	EA Price ₩
G11-220-535	202238-B	88 cc	H <sub>2</sub> O: 6.3 g (Max. 11.0 g)	1/8" Brass Compression	247,800
G11-220-541	202239-B			1/4" Brass Compression	247,800
G11-220-502	202238-4-PC			1/4" Push-to-connect	247,800
G11-220-506	202238-6-PC			6mm Push-to-connect	247,800
G11-220-510	202238L-B	116 cc	H <sub>2</sub> O: 8.3 g (Max. 14.4 g)	1/8" Brass Compression	273,100
G11-220-514	202239L-B			1/4" Brass Compression	273,100
G11-220-518	202238L-4-PC			1/4" Push-to-connect	273,100
G11-220-522	202238L-6-PC			6mm Push-to-connect	273,100
G11-220-526	202238D-B	240 cc	H <sub>2</sub> O: 17.2 g (Max. 29.9 g)	1/8" Brass Compression	273,100
G11-220-530	202239D-B			1/4" Brass Compression	273,100
G11-220-534	202238D-4-PC			1/4" Push-to-connect	273,100
G11-220-538	202238D-6-PC			6mm Push-to-connect	273,100
G11-220-547	202238XL-B	400 cc	H <sub>2</sub> O: 28.6 g (Max. 49.8 g)	1/8" Brass Compression	330,500
G11-220-563	202239XL-B			1/4" Brass Compression	330,500
G11-220-550	202238XL-4-PC			1/4" Push-to-connect	330,500
G11-220-554	202238XL-6-PC			6mm Push-to-connect	330,500
Capacity	88 cc	116 cc	240 cc	400 cc	
Nominal Flow Rate	0.68 SLPM	0.9 SLPM	1.85 SLPM	3.0 SLPM	
Max. Flow Rate	6.8 SLPM	9.0 SLPM	18.5 SLPM	30.0 SLPM	



### SPure™ HC Purifier, 플라스틱 케이스 탄화수소 트랩

Cat. No.	Model	용량	Trap 용량	Fitting	EA Price ₩
G11-220-600	202326-B	88 cc	HC: 7.2 g (Max. 24.2 g)	1/8" Brass	247,800
G11-220-603	202327-B			1/4" Brass	247,800
G11-220-608	202326-4-PC			1/4" Push	273,100
G11-220-611	202326-6-PC			6mm Push	273,100
G11-220-616	202326L-B	116 cc	HC: 9.6 g (Max. 31.9 g)	1/8" Brass	273,100
G11-220-621	202327L-B			1/4" Brass	273,100
G11-220-624	202326L-4-PC			1/4" Push	273,100
G11-220-628	202326L-6-PC			6mm Push	273,100
G11-220-632	202326D-B	240 cc	HC: 19.8 g (Max. 66.0 g)	1/8" Brass	273,100
G11-220-636	202327D-B			1/4" Brass	273,100
G11-220-641	202326D-4-PC			1/4" Push	273,100
G11-220-644	202326D-6-PC			6mm Push	273,100
G11-220-648	202326XL-B	400 cc	HC: 33.0 g (Max. 110.0 g)	1/8" Brass	330,500
G11-220-652	202327XL-B			1/4" Brass	330,500
G11-220-656	202326XL-4-PC			1/4" Push	330,500
G11-220-660	202326XL-6-PC			6mm Push	330,500
Capacity	88 cc	116 cc	240 cc	400 cc	
Nominal Flow Rate	0.68 SLPM	0.9 SLPM	1.85 SLPM	3.0 SLPM	
Max. Flow Rate	3.4 SLPM	4.5 SLPM	9.25 SLPM	15.0 SLPM	

### SPure™ H<sub>2</sub>O / HC Purifier, 수분 / 탄화수소 트랩 Indicator : Green to Blue (at 5% Relative Humidity)

Cat. No.	Model	용량	Trap 용량	Fitting	EA Price ₩
G11-220-664	202324-B	88 cc	H <sub>2</sub> O: 4.2 g (Max. 7.3 g) / HC: 2.9 g (Max. 9.7 g)	1/8" Brass	247,800
G11-220-668	202325-B			1/4" Brass	247,800
G11-220-671	202324-4-PC			1/4" Push	247,800
G11-220-674	202324-6-PC			6mm Push	247,800
G11-220-679	202324L-B	116 cc	H <sub>2</sub> O: 5.5 g (Max. 9.6 g) / HC: 3.8 g (Max. 12.8 g)	1/8" Brass	273,100
G11-220-683	202325L-B			1/4" Brass	273,100
G11-220-687	202324L-4-PC			1/4" Push	273,100
G11-220-672	202324L-6-PC			6mm Push	273,100
G11-220-676	202324D-B	240 cc	H <sub>2</sub> O: 11.5 g (Max. 19.9 g) / HC: 7.9 g (Max. 26.5 g)	1/8" Brass	273,100
G11-220-680	202325D-B			1/4" Brass	273,100
G11-220-684	202324D-4-PC			1/4" Push	273,100
G11-220-688	202324D-6-PC			6mm Push	273,100
G11-220-692	202324XL-B	400 cc	H <sub>2</sub> O: 19.1 g (Max. 33.2 g) / HC: 13.2 g (Max. 44.1 g)	1/8" Brass	330,500
G11-220-696	202325XL-B			1/4" Brass	330,500
G11-220-705	202324XL-4-PC			1/4" Push	330,500
G11-220-710	202324XL-6-PC			6mm Push	330,500
Capacity	88 cc	116 cc	240 cc	400 cc	
Nominal Flow Rate	0.68 SLPM	0.9 SLPM	1.85 SLPM	3.0 SLPM	
Max. Flow Rate	3.4 SLPM	4.5 SLPM	9.25 SLPM	15.0 SLPM	