

ZPure™ Inline Gas Filter, Stainless Steel

고효율 가스 정제관

- 고성능 Activated Absorbent와 Long Purifier Life
- 다양한 필터 용량 : 130cc, 475cc, 500cc, 750cc, 3.8L
- All Stainless Steel Body 와 Stainless Steel Compression Fitting
- Optional Fitting Configuration : 요청에의해 공급 가능
 - Compression Fitting : Max. 1000psi, Brass (1/4" or 1/8")
 - Quick Connect Fitting : Max.200psi, SS or Brass (1/4" or 1/8")
 - Push-to-Connect Fitting : Max.150psi, Brass (1/4" or 6mm)

• Chromres

Capacity	130cc	475cc	500cc	750cc	3.8L
Size, mm	φ32 x h280	φ38 x h570	φ50 x h350	φ50 x h500	φ100 x h640
Max. Pressure	68.9 bar / 1000 psi		13.8 bar / 200 psi		150 psi Air, N2 80 psi H2, O2

ZPure™ H₂O Purifier, 고효율 수분 트랩 ①

Cat. No.	Model	용량	Trap 용량 / Efficiency	Fitting	EA Price ₩
G11-220-755	202230-SS	130cc	H ₂ O: 12 g (Max. 22 g) / < 20 ppb	1/8"	339,200
G11-220-761	202231-SS			1/4"	339,200
G11-220-715	202230L-SS	475cc	H ₂ O: 45 g (Max. 79 g) / < 20 ppb	1/8"	500,600
G11-220-720	202231L-SS			1/4"	500,600
G11-220-725	202230D-SS	500cc	H ₂ O: 48 g (Max. 83 g) / < 20 ppb	1/8"	500,600
G11-220-730	202231D-SS			1/4"	500,600
G11-220-735	202230XL-SS	750cc	H ₂ O: 72 g (Max. 124 g) / < 20 ppb	1/8"	558,800
G11-220-740	202231XL-SS			1/4"	558,800
G11-220-745	202231M-SS	3.8L	H ₂ O: 369 g (Max. 640 g) / < 5 ppb	1/4"	5,914,000
G11-220-750	202230M-V			1/2" VCR	6,248,700

Capacity	130cc	475cc	500cc	750cc	3.8L
Nominal Flow Rate	1 SLPM	3.7 SLPM	3.8 SLPM	5.8 SLPM	30 SLPM
Max. Flow Rate	10 SLPM	36 SLPM	38 SLPM	57 SLPM	200 SLPM

Zpure O₂ Purifier, 고효율 산소 트랩 ②

Cat. No.	Model	용량	Trap 용량 / Efficiency	Fitting	EA Price ₩
G11-220-711	202216-SS	130cc	O ₂ 2500 cc / < 5 ppb	1/8"	681,900
G11-220-717	202217-SS			1/4"	681,900

- Nominal Flow Rate : 1 SLPM, Max. Flow Rate : 1 SLPM

ZPure™ HC Purifier, 고효율 탄화수소 트랩 ③

Cat. No.	Model	용량	Trap 용량 / Efficiency	Fitting	EA Price ₩
G11-220-701	202232-SS	130cc	HC: 11 g (Max. 36 g) / < 5 ppb	1/8"	339,200
G11-220-708	202233-SS			1/4"	339,200
G11-220-721	202232L-SS	475cc	HC: 40 g (Max. 131 g) / < 5 ppb	1/8"	500,600
G11-220-728	202233L-SS			1/4"	500,600
G11-220-731	202232D-SS	500cc	HC: 42 g (Max. 137 g) / < 5 ppb	1/8"	500,600
G11-220-737	202233D-SS			1/4"	500,600
G11-220-741	202232XL-SS	750cc	HC: 63 g (Max. 206 g) / < 5 ppb	1/8"	558,800
G11-220-748	202233XL-SS			1/4"	558,800
G11-220-751	202233M-SS	3.8L	HC: 327 g (Max. 1060 g) / < 5 ppb	1/4"	5,914,000
G11-220-758	202232M-V			1/2" VCR	6,248,700

Capacity	130cc	475cc	500cc	750cc	3.8L
Nominal Flow Rate	1 SLPM	3.7 SLPM	3.8 SLPM	5.8 SLPM	30 SLPM
Max. Flow Rate	10 SLPM	36 SLPM	38 SLPM	57 SLPM	200 SLPM

ZPure™ CO₂ Purifier, 고효율 이산화탄소 트랩 ④

Cat. No.	Model	용량	Trap 용량 / Efficiency	Fitting	EA Price ₩
G11-220-703	202212-SS	130cc	CO ₂ : 13 g (Max. 19 g) / < 5 ppb	1/8"	634,600
G11-220-709	202213-SS			1/4"	634,600
G11-220-713	202212L-SS	475cc	CO ₂ : 47 g (Max. 69 g) / < 5 ppbb	1/8"	1,113,600
G11-220-719	202213L-SS			1/4"	1,113,600
G11-220-723	202212D-SS	500cc	CO ₂ : 50 g (Max. 73 g) / < 5 ppb	1/8"	1,113,600
G11-220-729	202213D-SS			1/4"	1,113,600
G11-220-733	202212XL-SS	750cc	CO ₂ : 60 g (Max. 110 g) / < 5 ppb	1/8"	1,250,800
G11-220-739	202213XL-SS			1/4"	1,250,800
G11-220-749	202212M-V	3.8L	CO ₂ : 300 g (Max. 564 g) / < 5 ppb	1/2" VCR	11,330,000

Capacity	130cc	475cc	500cc	750cc	3.8L
Nominal Flow Rate	650 cc/min	2.5 SLPM	2.5 SLPM	4 SLPM	< 20 SLPM



다음 페이지에서 계속...