


ZPure™ Glass Indicating Purifier, Visual Indicator 에서 계속...
ZPure™ Glass H₂O / HC Purifier & Indicator, 멀티 가스용 수분 지시관 Green to Blue (Moisture)

Cat. No.	Model	용량	Trap 용량	Efficiency	Fitting	EA Price ₩
G11-220-060	202256-B	46 cc	H ₂ O : 3 g HC : 2 g	H ₂ O : < 20 ppb HC : < 5 ppb	1/8" Brass	314,300
G11-220-068	202257-B				1/4" Brass	314,300
G11-220-070	202256-SS				1/8" SS	415,200
G11-220-078	202257-SS				1/4" SS	415,200

- Dry Air or Inert Gas로부터 Moisture and Hydrocarbon 의 제거
- Gases to be Purified : Inert Carrier Gas, He, N₂, H₂, Ar, Air


Wall Mounting for CRS Filter Holder

Cat. No.	Model	Description	EA Price ₩
G11-220-100	205705	Mounting Kit for 1 Purifier, 2 x Clip	31,800

- 130cc SS ZPure Filter, ZPure Glass Indicating Filter, and 88cc or 116cc SPure Plastic Filter용
- 구성 : 2 Clips, Aassociated Hardware for Mounting (for 1 Gas Filter)

Type	Polygas I	Polygas II	O ₂	H ₂ O HSMI	H ₂ O	H ₂ O/HC
Normal Flow Rate	200 ml/miin					
Max. Flow Rate	400 ml/min		2 SLPM			
Max. Pressure	68.9 bar / 1000 psi					
Size	ø38 x h263 mm					


Advanced Filter System, AFS Purifier
가스 정제 장치


Cat. No.	Model	Trap 용량 (Efficiency)	Fitting	EA Price ₩
G11-220-007	202900	O ₂ : 850cc (< 5ppb), H ₂ O : 12g (<20ppb), HC : 8g (<5ppb)	SS 1/8"	2,672,500
G11-220-010	202902		SS 1/4"	2,672,500
G11-220-014	202903	H ₂ O : 8.8g (<20ppb), HC : 13g (<5ppb)	SS 1/8"	2,672,500
G11-220-018	202904		SS 1/4"	2,672,500
G11-220-022	202905	H ₂ O : 20g (<20ppb)	SS 1/8"	2,672,500
G11-220-026	202906		SS 1/4"	2,672,500

- 하나의 카트리지로 대용량 및 여러가스의 효과적인 가스 정제가능
- GC, GC/MS Carrier Gas, ECD Make up Gas, ICP에 적용 가능
- 구성 : SS Manifold, 대용량 카트리지, SS Fitting, Bracket for Wall or Bench
- Metal과 Glass로만 구성되며 Visible Indicator 장착
- GC와 1년 또는 그 이상의 사용 가능
- 카트리지 교환이 간단함 (교환시 Check Valve 작동으로 오염되지 않음)
- Oxygen : Green to Grey, Moisture : Blue to Brown
- **Chromres**

Normal Flow Rate	0.5 SLPM (Max. Recommended Flow Rate : 2 SLPM)
Max. Pressure	13.8 bar / 200 pai
Size	300 x 50 x 50 mm

Replacement Cartridge only

Cat. No.	Model	Description	EA Price ₩
G11-220-911	202910	for O ₂ , H ₂ O, HC	1,702,400
G11-220-921	202950	for H ₂ O, HC	1,106,600
G11-220-931	202955	for H ₂ O	1,106,600


Refill Kit for "Chromres" Gas Purifier
가스 정제관용 충전제

Cat. No.	Model	Description	EA Price ₩
G11-220-511	202298	CoFree(TM) Indicator, 350cc / Molecular Sieve 5A, 120cc	158,900
G11-220-411	202247	CoFree(TM) Indicator 500cc, Moisture 제거용	158,900
G11-220-311	202245	Molecular Sieve 5A 500cc, 불순물 제거용	158,900
G11-220-211	202243	Processed Activated Carbon 500cc	158,900

- Plastic Bottle로 포장
- **Chromres**

