


Blank Round Bottom Flask 초자 가공용 환저 플라스크, LukeGL®

Cat. No.	Model	용량	Size	EA Price ₩
J02-503-321	7000050	50 ml	∅51×h105 mm	2,400
J02-503-325	7000100	100 ml	∅64×h110 mm	2,500
J02-503-329	7000250	250 ml	∅85×h145 mm	3,900
J02-503-333	7000500	500 ml	∅105×h168 mm	5,700
J02-503-337	7001000	1 L	∅131×h210 mm	9,000
J02-503-341	7002000	2 L	∅166×h260 mm	15,900
J02-503-343	7003000B	3 L	∅185×h260 mm	48,600


- Borosilicate Glass 3.3
- DIN12347, ISO1773 규격에 따름
- Blank Type으로 초자 가공용으로 적합
- LukeGL®



Evaporating Flask, Blank 초자 가공용 농축 플라스크, LukeGL®

Cat. No.	Model	용량	Size	EA Price ₩
J02-503-555	7100050	50 ml	∅51×h105 mm	8,900
J02-503-559	7100100	100 ml	∅64×h110 mm	10,700
J02-503-563	7100250	250 ml	∅85×h145 mm	14,500
J02-503-567	7100500	500 ml	∅105×h168 mm	18,800
J02-503-571	7101000	1 L	∅131×h210 mm	38,000
J02-503-575	7102000	2 L	∅166×h260 mm	50,100

- Borosilicate Glass 3.3
- Blank Type으로 초자 가공용으로 적합
- DIN12347, ISO1773 규격에 따름
- With Beaded Rim
- LukeGL®



Round Bottom Flask, Simax® 환저 플라스크

Cat. No.	Model	용량	Size	Neck Dia.	EA Price ₩
J02-985-003	205050	50 ml	∅51×h105 mm	22 mm	
J02-985-006	205100	100 ml	∅64×h115 mm	22 mm	
J02-985-009	205250	250 ml	∅85×h145 mm	34 mm	
J02-985-012	205500	500 ml	∅105×h175 mm	34 mm	
J02-985-015	205940	1 L	∅131×h210 mm	42 mm	
J02-985-018	205950	2 L	∅166×h260 mm	42 mm	
J02-985-021	205954	4 L	∅207×h315 mm	50 mm	
J02-985-024	205956	5 L	∅222×h330 mm	50 mm	
J02-985-027	205958	6 L	∅236×h355 mm	65 mm	
J02-985-030	205966	10 L	∅279×h420 mm	65 mm	
J02-985-033	205968	20 L	∅345×h500 mm	80 mm	

- Simax® Borosilicate Glass 3.3
- Neck 부위가 약간 사선으로 되어있음
- KAVALIERRGLASS

