


JY Spherical Joint, O-ring Type 오링형 볼 조인트

- Grease-free Type의 볼 조인트
- 고 진공용으로 사용하기에 알맞음
- J Young Brand

Greaseless Ball Joint, Without O-ring ①

Cat. No.	Model	Bore	Shark od	EA Price ₩
J02-179-103	GBM/12/2	2 mm	7	44,100
J02-179-105	GBM/12/5	5 mm	8	32,000
J02-179-107	GBM/18/7	7 mm	11	34,800
J02-179-109	GBM/18/9	9 mm	12	34,800
J02-179-111	GBM/28/12	12 mm	17	44,600
J02-179-113	GBM/28/15	15 mm	19	44,600
J02-179-115	GBM/35/20	20 mm	24	51,600
J02-179-117	GBM/40/25	25 mm	28	67,000
J02-179-119	GBM/51/35	35 mm	38	86,900



Greaseless Cup Joint ②

Cat. No.	Model	Bore	Shark od	EA Price ₩
J02-179-123	GBF/12/2	2 mm	7	44,100
J02-179-125	GBF/12/5	5 mm	8	32,000
J02-179-127	GBF/18/7	7 mm	11	34,800
J02-179-129	GBF/18/9	9 mm	12	34,800
J02-179-131	GBF/28/12	12 mm	17	44,600
J02-179-133	GBF/28/15	15 mm	19	44,600
J02-179-135	GBF/35/20	20 mm	24	51,600
J02-179-137	GBF/40/25	25 mm	28	67,000
J02-179-139	GBF/51/35	35 mm	38	86,900



JYO-ring for Spherical Joint

Viton O-ring

Cat. No.	Model	EA Price ₩
J02-179-143	O12/V	1,600
J02-179-145	O18/V	1,600
J02-179-147	O28/V	2,000
J02-179-149	O35/V	2,400
J02-179-151	O40/V	2,700

FEP Covered Viton O-ring ③

Cat. No.	Model	EA Price ₩
J02-179-155	O12/P	57,000
J02-179-157	O18/P	86,500
J02-179-159	O28/P	129,800
J02-179-161	O35/P	188,900
J02-179-163	O40/P	183,700

적용 규격
12/2 & 12/5
18/7 & 18/9
28/12 & 28/15
35/20
40/25

JY Clip for Spherical Joint


Thumb Screw Lock Clip ④

Cat. No.	Model	EA Price ₩
J02-179-183	TJC/13	69,800
J02-179-185	TJC/19	80,500
J02-179-187	TJC/29	80,500
J02-179-189	TJC/35	111,200
J02-179-191	TJC/41	151,000

적용 규격
12/2 & 12/5
18/7 & 18/9
28/12 & 28/15
35/20
40/25

Gleazeless Conical Joint, O-ring Type 고진공용 유리 조인트, Gleaseless

- Grease-free Type의 유리 텀퍼 조인트
- Glass Cone은 Viton O-ring이 부착
- Clear Glass Socket
- Vacuum Application : 10⁻⁶ Torr • J Young Brand



Cone

Cat. No.	Model	EA Price ₩
J02-179-063	ORB/10	32,000
J02-179-065	ORB/14	34,100
J02-179-067	ORB/19	36,100
J02-179-069	ORB/24	54,200
J02-179-071	ORB/29	54,600
J02-179-073	ORB/34	60,700

Socket with Viton O-ring

Cat. No.	Model	EA Price ₩
J02-179-083	ORS/10	31,100
J02-179-085	ORS/14	33,300
J02-179-087	ORS/19	35,100
J02-179-089	ORS/24	43,400
J02-179-091	ORS/29	47,300
J02-179-093	ORS/34	59,500

적용 규격
B10
B14
B19
B24
B29
B34