

General Purpose Needle Valve / Stopcock

테프론 니들 밸브

- Self Retaining Key, No Ground Surface
- 적정 사용 온도 : -18 ~ 60°C
- Never Leak & Seize
- J Young Brand

Straight Type, 일자형 니들 밸브 ①



Cat. No.	Model	Side Arm id	Side Arm od	EA Price ₩
S11-179-301	GP/5	5 mm	8 mm	89,200
S11-179-303	GP/8	8 mm	11 mm	89,300

Right Angle Type, 90° 형 니들 밸브 ②

Cat. No.	Model	Side Arm id	Side Arm od	EA Price ₩
S11-179-309	GP/5/RA	5 mm	8 mm	73,400
S11-179-311	GP/8/RA	8 mm	11 mm	73,500

Double Pattern, 더블형 니들 밸브 ③

Cat. No.	Model	Side Arm id	Side Arm od	EA Price ₩
S11-179-317	GP/D/5	5 mm	8 mm	135,500
S11-179-319	GP/D/10	2.7 mm	7 mm	188,400




PTFE High Vacuum Screw Cap Joint

고진공용 스크류 캡 연결관

Cat. No.	Model	Type	Thread	Tubing id	Tubing od	PTFE Insert Bore	EA Price ₩
J02-179-203	DSY/0/1	일자형	13 mm	5 mm	8 mm	8 mm	81,700
J02-179-205	DSY/0/2		19 mm	8 mm	11 mm	11 mm	81,900
J02-179-207	DSY/00/1		13 mm	5 mm	8 mm	8 mm	102,300
J02-179-209	DSY/00/2		19 mm	8 mm	11 mm	11 mm	102,800
J02-179-211	DSY/1/1	L자형	13 mm	5 mm	8 mm	8 mm	66,700
J02-179-213	DSY/1/2		19 mm	8 mm	11 mm	11 mm	61,000
J02-179-215	DSY/2/1	T자형	13 mm	5 mm	8 mm	8 mm	128,100
J02-179-217	DSY/2/2		19 mm	8 mm	11 mm	11 mm	129,300
J02-179-219	DSY/3/1	십자형	13 mm	5 mm	8 mm	8 mm	159,200
J02-179-221	DSY/3/2		19 mm	8 mm	11 mm	11 mm	160,000
J02-179-223	DSY/4/1	십자형	13 mm	5 mm	8 mm	8 mm	211,000
J02-179-225	DSY/4/2		19 mm	8 mm	11 mm	11 mm	211,800

- J Young Stopcock (POR/5 & POR/8)과 직접 연결 가능함
- Screw Cap을 풀어 연결되어진 Tubing을 조절 할 수 있음
- Glass Insert Tubing, PTFE Insert & Screw Cap 포함
- J Young Brand



High Vacuum Needle Valve

고진공 니들 밸브

Cat. No.	Model	Type	Bore	EA Price ₩
S11-132-011	963-01	90° 형 ①	0 ~ 4 mm	110,000
S11-132-018	965-01	좌우 일자형 ②		122,400
S11-132-025	967-01	상하 일자형 ③		125,700

- 일반형 보다 정밀하게 유량을 조절할 수 있음
- Bore Size : 0 ~ 4 mm
- Chemglass

General Purpose Valve for Liquid

액체용 밸브

Cat. No.	Model	Type	Bore	EA Price ₩
S11-132-121	587-01	90° 형 ①	2 mm	92,700
S11-132-131	591-01	좌우 일자형 ②		92,200
S11-132-141	595-01	상하 일자형 ③		92,200

- O-ring이 없는 Teflon Plug와 유리로 구성
- Semi-needle Valve와 같은 기능을 함
- Chemglass

