

High Vacuum Valve, 1-arm 고진공 밸브

- Vacuum Application : Down to 10^{-7} torr
- Single sidearm sealed at a 90° angle to the barrel
- 일체형 Teflon Tip O-ring 및 Teflon Coating Viton O-ring을 통한 Sealing
- 유체와 닿는 부분은 모두 유리나 Teflon Part • Chemglass

Cat. No.	Model	Bore	EA Price ₩
S11-132-211	960-01	0 ~ 4 mm	85,900
S11-132-216	960-02	0 ~ 8 mm	132,000
S11-132-221	960-03	0 ~ 12 mm	136,500
S11-132-226	960-04	0 ~ 20 mm	457,300



High Vacuum Valve 고진공 밸브

- Vacuum Application : Down to 10^{-7} torr
- 일체형 Teflon Tip O-ring 및 Teflon Coating Viton O-ring을 통한 Sealing
- 유체와 닿는 부분은 모두 유리나 Teflon Part • Chemglass

90° 형 ①

Cat. No.	Model	Bore	EA Price ₩
S11-132-621	951-01	0 ~ 4 mm	102,500
S11-132-628	951-02	0 ~ 8 mm	122,000
S11-132-635	951-03	0 ~ 12 mm	141,900

좌우 일자형 ②

Cat. No.	Model	Bore	EA Price ₩
S11-132-792	953-01	0 ~ 4 mm	115,500
S11-132-799	953-02	0 ~ 8 mm	168,200
S11-132-806	953-03	0 ~ 12 mm	185,700

상하 일자형 ③

Cat. No.	Model	Bore	EA Price ₩
S11-132-459	955-01	0 ~ 4 mm	117,200
S11-132-466	955-02	0 ~ 8 mm	181,300

3-way 형 ④

Cat. No.	Model	Bore	EA Price ₩
S11-132-888	957-01	0 ~ 4 mm	243,600
S11-132-895	957-02	0 ~ 8 mm	314,700



Pressure Relief Valve 압력 조절 밸브

Cat. No.	Model	Relief Pressure	EA Price ₩
S11-132-185	999-01	5 ~ 50 psi	136,300

- Stainless Steel Spring의 Tension에 따른 압력 조절형
- Automatic, Adjustable Pressure Relief Valve
- Chemglass



Inlet Valve 진공 투입 및 벤트 밸브

Cat. No.	Model	Bore	EA Price ₩
S11-132-199	590-01	4 mm	80,100

- 일자형 Stopcock의 대체용 Valve
- Teflon Coated Viton O-ring
- Chemglass

