

Pharma Plus-A Glove Box Glove /
Ethylene-Propylene-Diene-Rubber (EPDM)
EPDM 글러브 박스 글러브

Cat. No.	Model	Hand size	Port sizes	EA Price ₩
G04-508-249	16BL05 180	L (9~10)	7 (7.5"~ 8.5")	1,252,200
G04-508-251	16BL05 230		9 (9.4" ~ 11")	1,252,200
G04-508-252	16BL05 250		10 (10.4" ~ 12")	1,287,200

- Thickness : 20mil (0.5mm), Length : 800mm
- FDA 적합 재료로 제작 • 자외선에 강함
- 고압멸균 후에도 끈적이지 않고 침습적이지 않음
- 과산화수소 및 일반적인 소독 물질에 대한 내화학성이 강함
- Halogen Free (소각 폐기 가능) • 제약, 식품 산업 등에서 사용 가능
- Black Liner로 인해 외부 손상을 쉽게 확인 가능

Pharma-A Glove Box Glove
Ethylene-Propylene-Diene-Rubber (EPDM)
EPDM 글러브 박스 글러브

Cat. No.	Model	Hand size	Port sizes	EA Price ₩
G04-508-199	17BL04 180	L (9~10)	7 (7.5"~ 8.5")	944,300
G04-508-201	17BL04 230		9 (9.4" ~ 11")	944,300
G04-508-202	17BL04 250		10 (10.4" ~ 12")	979,400

- Thickness : 16mil (0.4mm), Length : 800mm
- FDA 적합 재료로 제작 • 편안한 착용감 • Halogen Free (소각 폐기 가능)
- 고압멸균 후에도 끈적이지 않고 침습적이지 않음
- 과산화수소 및 일반적인 소독 물질에 대한 내화학성이 강함
- 자외선에 강함 • 106 Ω 이상의 저항으로 접지시 정전기가 없음
- 영구변형이 적음 • 내부압력 증가 시 빠르게 팽창해 빠른 누출 확인 가능
- 제약 및 생명공학 산업 등에서 사용 가능

Isflex Glove Box Glove / XSBR Elastomer
이소플렉스 글러브 박스 글러브

Cat. No.	Model	Thickness	Hand size	Port Sizes	EA Price ₩
G04-508-014	10BL05 180	20mil (0.5mm)	L (9~10)	7.5-8.5"	746,900
G04-508-016	10BL05 230			9.4~11"	746,900
G04-508-017	10BL05 250			10.4 ~ 12"	781,700

- 과산화수소 / 아이소프로필 알코올에 사용 가능
- FDA 적합 재료로 제작 • Latex Free
- 착용감이 좋고 민감한 작업에 적합 • 오존 및 마모에 내구성이 좋음
- 제약, 의료 및 식품 산업에 적합 • 제약 산업의 격리환경에 적합



Jugitec® 글러브 박스용 글러브 상세 스펙

Specifcaions (선택사항)

Polymer Type	Bromobutyl-Rubber (BIIR)	Hypalon (CSM)	Ethylene-Propylene-Diene-Rubber (EPDM) Pharma Plus-A	Ethylene-Propylene-Diene-Rubber (EPDM) Pharma-A	XSBR-Elastomer
Color	Black	White	White (Black Liner)	Black	White
Hand Size	L (9~10) / XL (11)	L (9~10) / XL (11)	L (9~10) / XL (11)	L (9~10) / XL (11)	L (9~10) / XL (11)
Length	800 mm / 920 mm	800 mm / 920 mm	800 mm / 920 mm	800 mm / 920 mm	800 mm / 920 mm
Hand types	fully anatomical / ambidextrous / tactile	fully anatomical / ambidextrous / tactile	fully anatomical / ambidextrous / tactile	fully anatomical / ambidextrous / tactile	fully anatomical / ambidextrous / tactile - available only in size L
Portsize	5.7" (5.9~6.73") 7" (7.5~8.5") 9" (9.4~11") 10" (10.4 ~ 12")	7" (7.5~8.5") 9" (9.4~11") 10" (10.4 ~ 12")	7" (7.5~8.5") 9" (9.4~11") 10" (10.4 ~ 12")	5.7" (5.9~6.73") 7" (7.5~8.5") 9" (9.4~11") 10" (10.4 ~ 12")	7" (7.5~8.5") 9" (9.4~11") 10" (10.4 ~ 12")
Thickness	16mil (0.4mm) 24mil (0.6mm) 31mil (0.8mm)	16mil (0.4mm) 24mil (0.6mm) 31mil (0.8mm)	20mil (0.5mm)	16mil (0.4mm) 24mil (0.6mm)	20mil (0.5mm)
Temperature Range	-40° C ~ 90° C	-20° C ~ 120° C	-20° C ~ 130° C Autoclave	-20° C ~ 130° C Autoclave	-20° C ~ 80° C