


## Compact Glove Box for Balance 저울용 글러브 박스, for Analytical & Micro Balance

배송비 별도

Cat. No.	Model	Size	EA Price ₩
G04-566-191	830-ABC	w1100×d610×h790 mm	17,333,900
G04-566-201	830-ABB	w1250×d760×h790 mm	19,084,900
G04-566-211	830-ABD	w1600×d790×h860 mm	20,003,000
G04-566-221	830-ABE	w1930×d812×h890 mm	21,553,200

- 외부와 차단된 환경에서(위험물, 수분에 민감한 시료 등) 계량을 위한 Glove Box (분석 저울, 마이크로 저울 등)
- 수용액의 수분량, 접착제, 곡물류, 독성 화학 물질, 의약품, 종이, 플라스틱, 소량의 방사선 활성 물질과 유해물질의 무게 계량에 적합
- 투명 아크릴 구조 (10mm 두께)
- 석정반 위에서 사용하는 것이 좋음
- 전자 저울의 이동을 위해 상부 탈착 가능
- 클램프의 마모를 최소화하기 위해 클램프 조정 가능
- HEPA Filter와 Carbon Filter가 결합된 830 Filter 장착
- 4개의 Key Cock Valve
- 백색의 Hypalon 글러브 사용
- UL, CSA, and CE approved
- PLAS LABS



Cat. No.	G04-566-191	G04-566-201	G04-566-211	G04-566-221	
Model	830-ABC	830-ABB	830-ABD	830-ABE	
Size	Internal	w710×d580×h740 mm	w890×d740×h760 mm	w1220×d740×h810 mm	w1524×d740×h810 mm
	External	w1100×d610×h790 mm	w1250×d760×h790 mm	w1600×d790×h860 mm	w1930×d812×h890 mm
	Top Opening	w432×d585 mm	w610×d735 mm	w610×d1065 mm	w610×d1370 mm
	Access Door Opening (Option)	w400×d560 mm	w546×d560 mm	w546×d610 mm	w546×d610 mm
Approx. Volume	306 L	498 L	730 L	912 L	
Approx. Weight	104 kg	125 kg	141 kg	182 kg	
Port	2ea	2ea	3ea	4ea	

## Accessory for Glove Box 글러브 박스용 악세서리

Cat. No.	Model	Description	EA Price ₩
G04-566-513	800-ONEG	Oxygen Reducing Accessory with Pump	3,282,000
G04-566-517	800-DT	Drying Train Package with Pump	4,200,100
G04-566-521	800-AS/SPI	Anti-Static Ionizer	1,838,900
G04-566-525	800-T	"T" Transfer Chamber	8,414,100
G04-566-529	800-PLUG	Glove Port Plug (2ea)	915,300

• PLAS LABS



- Glove Box용 Glove는 566 ~ 567 Page 참조

