

Standard Square Grid 스퀘어 그리드



• Gilder

재질 : Copper

Cat. No.	Model	Mesh	Pitch	Width		PK Price ₩	
				Bar	Hole		
G05-167-004	G50-C	50	500 um	80	420	100ea	33,500
G05-167-008	G100-C	100	250 um	45	205	100ea	33,500
G05-167-012	G150-C	150	165 um	40	125	100ea	36,300
G05-167-016	G200-C	200	125 um	35	90	100ea	36,300
G05-167-040	G300-C	300	83 um	25	58	100ea	36,300
G05-167-044	G400-C	400	62 um	25	37	100ea	39,200
G05-167-048	G600-C	600	42 um	16 ~ 10	30	100ea	62,200

재질 : Nickel

Cat. No.	Model	Mesh	Pitch	Width		PK Price ₩	
				Bar	Hole		
G05-167-304	G50-N	50	500 um	80	420	100ea	46,600
G05-167-308	G100-N	100	250 um	45	205	100ea	46,600
G05-167-312	G150-N	150	165 um	40	125	100ea	50,600
G05-167-316	G200-N	200	125 um	35	90	100ea	50,600
G05-167-340	G300-N	300	83 um	25	58	100ea	50,600
G05-167-344	G400-N	400	62 um	25	37	100ea	51,000
G05-167-348	G600-N	600	42 um	16 ~ 10	30	100ea	77,300



Fine Square Mesh Grid HS 스퀘어 그리드

• Standard Grid 보다 매우 얇은 Bar를 가지고 있어 Open Area가 약 40% 더 높음 (Ultra-high Transmission)

• Gilder

재질 : Copper

Cat. No.	Model	Mesh	Pitch	Width		PK Price ₩	
				Bar	Hole		
G05-167-052	G200HS-C	200	125 um	12	113	100ea	64,300
G05-167-056	G300HS-C	300	83 um	10	73	100ea	74,500
G05-167-060	G400HS-C	400	62 um	8	54	100ea	97,100
G05-167-064	G1000HS-C	1000	25 um	6	19	25ea	68,200

재질 : Nickel

Cat. No.	Model	Mesh	Pitch	Width		PK Price ₩	
				Bar	Hole		
G05-167-352	G200HS-N	200	125 um	12	113	100ea	81,300
G05-167-356	G300HS-N	300	83 um	10	73	100ea	95,800
G05-167-360	G400HS-N	400	62 um	8	54	100ea	122,700
G05-167-364	G1000HS-N	1000	25 um	6	19	25ea	101,600



Slots Grid 슬롯 그리드

• Gilder

재질 : Copper

Cat. No.	Model	Slot Width	PK Price ₩
G05-167-172	GS2X0.5-C	2000 / 500 um	100ea 47,500
G05-167-176	GS2X1-C	2000 / 1000 um	100ea 47,500

재질 : Nickel

Cat. No.	Model	Slot Width	PK Price ₩
G05-167-472	GS2X0.5-N	2000 / 500 um	100ea 67,200
G05-167-476	GS2X1-N	2000 / 1000 um	100ea 67,200

