

Hexagonal Mesh Grid 헥사고널 그리드


- Standard Square Mesh와 같은 Pattern이나 Mesh의 모양이 육각으로 전자 현미경에서 좀 더 정확한 방향성을 제시
- Gilder

재질 : Copper

Cat. No.	Model	Mesh	Pitch	Width		PK Price ₩	
				Bar	Hole		
G05-167-112	G50HEX-C	50	500 um	70	430	100ea	33,500
G05-167-116	G75HEX-C	75	340 um	50	290	100ea	33,500
G05-167-120	G100HEX-C	100	250 um	35	215	100ea	33,500
G05-167-124	G150HEX-C	150	165 um	35	130	100ea	36,300
G05-167-128	G200HEX-C	200	125 um	25	100	100ea	36,300
G05-167-132	G300HEX-C	300	83 um	25	58	100ea	36,300
G05-167-136	G400HEX-C	400	62 um	25	37	100ea	39,200

재질 : Nickel

Cat. No.	Model	Mesh	Pitch	Width		PK Price ₩	
				Bar	Hole		
G05-167-412	G50HEX-N	50	500 um	70	430	100ea	46,600
G05-167-416	G75HEX-N	75	340 um	50	290	100ea	46,600
G05-167-420	G100HEX-N	100	250 um	35	215	100ea	46,600
G05-167-424	G150HEX-N	150	165 um	35	130	100ea	50,600
G05-167-428	G200HEX-N	200	125 um	25	100	100ea	50,600
G05-167-432	G300HEX-N	300	83 um	25	58	100ea	50,600
G05-167-436	G400HEX-N	400	62 um	25	37	100ea	51,000



Aperture Grid 싱글 홀 그리드


- Gilder

재질 : Copper

Cat. No.	Model	Hole Width	PK Price ₩
G05-167-150	GA100-C	100 um	100ea 53,500
G05-167-152	GA200-C	200 um	100ea 47,500
G05-167-156	GA300-C	300 um	100ea 47,500
G05-167-160	GA400-C	400 um	100ea 47,500
G05-167-164	GA500-C	500 um	100ea 47,500
G05-167-168	GA600-C	600 um	100ea 47,500

재질 : Nickel

Cat. No.	Model	Hole Width	PK Price ₩
G05-167-450	GA100-N	100 um	100ea 72,900
G05-167-452	GA200-N	200 um	100ea 67,200
G05-167-456	GA300-N	300 um	100ea 67,200
G05-167-460	GA400-N	400 um	100ea 67,200
G05-167-464	GA500-N	500 um	100ea 67,200
G05-167-468	GA600-N	600 um	100ea 67,200



Finder Grid 파인더 그리드, F4 Type

- Gilder

재질 : Copper

Cat. No.	Model	Mesh	Pitch	Width		PK Price ₩	
				Bar	Hole		
G05-167-218	G200F4-C	200/250	125/105	20/15	105/ 90	100ea	62,400

재질 : Nickel

Cat. No.	Model	Mesh	Pitch	Width		PK Price ₩	
				Bar	Hole		
G05-167-518	G200F4-N	200/250	125/105	20/15	105/ 90	100ea	95,200

