

## Finder Grid

## 파인더 그리드, SEMF1 Type

- Gilder

재질 : Copper

Cat. No.	Model	Mesh	Pitch	Width		PK Price W
				Bar	Hole	
G05-167-222	SEMF1-C	10	15	150	N/A	10ea 130,500

재질 : Nickel

Cat. No.	Model	Mesh	Pitch	Width		PK Price W
				Bar	Hole	
G05-167-522	SEMF1-N	10	15	150	N/A	10ea 154,100



## Finder Grid

## 파인더 그리드, SEMF2 Type


- Overall thickness : 50 um
- Large cells are numbered 1 to 57
- Gilder

재질 : Copper

Cat. No.	Model	Mesh	Pitch	Width		PK Price W
				Bar	Hole	
G05-167-226	SEMF2-C	10	500/ 1000	80/ 120	420/ 880	10ea 130,500

재질 : Nickel

Cat. No.	Model	Mesh	Pitch	Width		PK Price W
				Bar	Hole	
G05-167-526	SEMF2-N	10	500/ 1000	80/ 120	420/ 880	10ea 154,100

Optional Accessory for Grid  
그리드용 약세사리

## Grid Storage Box



Cat. No.	Model	Description	EA Price W
G05-167-911	SB50 ①	Store 50 Grids	21,300
G05-167-915	SB50BN	Store 50 Grids, with Unique Box #	24,000
G05-167-917	SB50RBN	Store 50 Grids, with Unique Box & Record Card	29,000
G05-167-920	SB100 ②	Store 100 Grids	24,000
G05-167-925	SB100BN ③	Store 100 Grids, with Unique Box #	27,100
G05-167-927	SB100RBN	Store 100 Grids, with Unique Box & Record Card	32,200

- SB50 Serise
  - Size : w75 × d65 × h6.5 mm
  - 무게 : 22 g
  - 재질 : ABS-PHAT(Base), CAB(Cover)
- SB100 Serise
  - Size : w75 × d65 × h6.5 mm
  - 무게 : 25 g
  - 재질 : ABS-PHAT(Base), CAB(Cover)
- Gilder

## Tweezer for Convenient Handling

Cat. No.	Model	Description	EA Price W
T14-33-088	RU-5-SA ④	Tweezer for handling grid	44,000

- Rubis



## Grid Handling용 Tweezer Set – High Precision Tweezer Page (1164 ~ 1171 Page 참조)

- Nano Tweezer
- General Type

