

Standard Weight Set, OIML 계속...

E2급 분동 세트 (In Plastic Case)

셋트구성 Value	포함수량 (Pcs)	Total 무게 (g)	in Plastic Case		EA Price ₩	DKD 인증서 포함 셋트		EA Price ₩
			Cat. No.	Model		Cat. No.	Model	
1 mg ~ 500 mg	12	1.11	W06-202-048	8.XFHP-710	768,800	W06-202-548	8.XFHP-711	1,294,800
1 mg ~ 5 g	16	6.11	W06-202-052	8.XFHP-730	1,182,300	W06-202-552	8.XFHP-731	1,604,700
1 mg ~ 50 g	20	111.11	W06-202-056	8.XFHP-760	1,711,700	W06-202-556	8.XFHP-761	2,239,800
1 mg ~ 100 g	21	211.11	W06-202-060	8.XFHP-770	1,858,400	W06-202-560	8.XFHP-771	2,430,800
1 mg ~ 200 g	23	611.11	W06-202-064	8.XFHP-780	2,220,200	W06-202-564	8.XFHP-781	2,881,200
1 mg ~ 500 g	24	1,111.11	W06-202-068	8.XFHP-790	3,019,400	W06-202-568	8.XFHP-791	3,213,600
1 mg ~ 1 kg	25	2,111.11	W06-202-072	8.XFHP-810	4,227,500	W06-202-572	8.XFHP-811	3,784,600
1 mg ~ 2 kg	27	6,111.11	W06-202-076	8.XFHP-820	5,127,500	W06-202-576	8.XFHP-821	5,025,300
1 mg ~ 5 kg	28	11,111.11	W06-202-080	8.XFHP-830	5,143,100	W06-202-580	8.XFHP-831	6,015,200
1 mg ~ 10 kg	29	21,111.11	W06-202-084	8.XFHP-840	6,495,900	W06-202-584	8.XFHP-841	7,429,300

- Wooden Box 분동 세트도 공급 가능함 (단, 1 mg ~ 10 kg Set는 제외)
- Tweezer 및 Dust Brush가 Plastic Case안에 기본 포함됨

“E2”구성 Component, 개별판매시 DKD 인증서 및 Plastic Case 기본 포함

Cat. No.	Model	Value	Tolerance 허용 오차	EA Price ₩	Cat. No.	Model	Value	Tolerance 허용 오차	EA Price ₩
W05-202-248	8.PAHP-011	1 mg	±0.006 mg	110,200	W05-202-748	8.MEHP-131	5 g	±0.050 mg	171,600
W05-202-252	8.PAHP-021	2 mg *	±0.006 mg	110,200	W05-202-752	8.MEHP-141	10 g	±0.060 mg	178,000
W05-202-256	8.PAHP-031	5 mg	±0.006 mg	110,200	W05-202-756	8.MEHP-151	20 g *	±0.080 mg	193,000
W05-202-260	8.PNHP-041	10 mg	±0.008 mg	110,200	W05-202-760	8.MEHP-161	50 g	±0.10 mg	214,300
W05-202-264	8.PNHP-051	20 mg *	±0.010 mg	110,200	W06-202-764	8.MEHP-171	100 g	±0.16 mg	314,000
W05-202-268	8.PNHP-061	50 mg	±0.012 mg	110,200	W06-202-768	8.MEHP-181	200 g *	±0.30 mg	352,800
W05-202-272	8.PNHP-071	100 mg	±0.016 mg	110,200	W06-202-772	8.MEHP-191	500 g	±0.80 mg	474,900
W05-202-276	8.PNHP-081	200 mg *	±0.020 mg	110,200	W06-202-776	8.MEHP-211	1 kg	±1.6 mg	746,300
W05-202-280	8.PNHP-091	500 mg	±0.025 mg	110,200	W06-202-780	8.MEHP-221	2 kg *	±3.0 mg	988,500
W06-202-284	8.MEHP-111	1 g	±0.030 mg	560,400	W06-202-784	8.MEHP-231	5 kg	±8.0 mg	1,298,700
W05-202-288	8.SEHP-121	2 g *	±0.040 mg	196,700	W06-202-788	8.MEHP-241	10 kg	±16 mg	1,831,700

- Value의 ‘*’ 의 표시 : 분동이 세트에 2개씩 포함됨
- 1mg ~ 500mg : 면분동

F1급 분동 세트 (In Plastic Box)

셋트구성 Value	포함수량 (Pcs)	Total 무게 (g)	in Plastic Box		EA Price ₩	DKD 인증서 포함 셋트		EA Price ₩
			Cat. No.	Model		Cat. No.	Model	
1 mg ~ 500 mg	12	1.11	W05-202-088	7.XFHP-710	523,500	W05-202-588	7.XFHP-711	647,600
1 mg ~ 5 g	16	6.11	W05-202-092	7.XFHP-730	829,800	W05-202-592	7.XFHP-731	1,015,700
1 mg ~ 50 g	20	111.11	W05-202-096	7.XFHP-760	1,217,700	W05-202-596	7.XFHP-761	1,465,400
1 mg ~ 100 g	21	211.11	W05-202-100	7.XFHP-770	1,328,800	W05-202-600	7.XFHP-771	1,599,200
1 mg ~ 200 g	23	611.11	W05-202-104	7.XFHP-780	1,598,600	W05-202-604	7.XFHP-781	1,919,500
1 mg ~ 500 g	24	1,111.11	W05-202-108	7.XFHP-790	1,851,000	W05-202-608	7.XFHP-791	2,206,500
1 mg ~ 1 kg	25	2,111.11	W05-202-112	7.XFHP-810	2,275,400	W05-202-612	7.XFHP-811	2,456,600
1 mg ~ 2 kg	27	6,111.11	W05-202-116	7.XFHP-820	3,294,700	W05-202-616	7.XFHP-821	3,549,300
1 mg ~ 5 kg	28	11,111.11	W05-202-120	7.XFHP-830	3,923,500	W05-202-620	7.XFHP-831	4,219,000
1 mg ~ 10 kg	29	21,111.11	W05-202-124	7.XFHP-840	5,080,400	W05-202-624	7.XFHP-841	5,842,400

- Wooden Box 분동 세트도 공급 가능함
- Tweezer 및 Dust Brush가 Plastic Box안에 기본 포함됨

“F1”구성 Component, 개별판매시 DKD 인증서 및 Plastic Case 기본 포함

Cat. No.	Model	Value	Tolerance 허용오차	EA Price ₩	Cat. No.	Model	Value	Tolerance 허용오차	EA Price ₩
W05-202-292	7.PAHP-011	1 mg	±0.020 mg	73,600	W05-202-793	7.MEHP-131	5 g	±0.16 mg	158,800
W05-202-296	7.PAHP-021	2 mg *	±0.020 mg	73,600	W05-202-797	7.MEHP-141	10 g	±0.20 mg	164,200
W05-202-300	7.PAHP-031	5 mg	±0.020 mg	73,600	W05-202-801	7.MEHP-151	20 g *	±0.25 mg	174,900
W05-202-304	7.PNHP-041	10 mg	±0.025 mg	73,600	W05-202-805	7.MEHP-161	50 g	±0.30 mg	196,500
W05-202-308	7.PNHP-051	20 mg *	±0.03 mg	73,600	W06-202-808	7.MEHP-171	100 g	±0.5 mg	244,300
W05-202-312	7.PNHP-061	50 mg	±0.04 mg	73,600	W06-202-812	7.MEHP-181	200 g *	±1.0 mg	271,400
W05-202-316	7.PNHP-071	100 mg	±0.05 mg	73,600	W06-202-816	7.MEHP-191	500 g	±2.5 mg	405,100
W05-202-320	7.PNHP-081	200 mg *	±0.06 mg	73,600	W06-202-820	7.MEHP-211	1 kg	±5 mg	600,900
W05-202-324	7.PNHP-091	500 mg	±0.08 mg	73,600	W06-202-824	7.MEHP-221	2 kg *	±10 mg	835,300
W05-202-329	7.MEHP-111	1 g	±0.10 mg	137,000	W06-202-828	7.MEHP-231	5 kg	±25 mg	971,000
W05-202-333	7.MEHP-121	2 g *	±0.12 mg	141,900	W06-202-832	7.MEHP-241	10 kg	±50 mg	1,571,900

- Value의 ‘*’ 의 표시 : 분동이 세트에 2개씩 포함됨
- 1mg ~ 500mg : 선분동