

PFA Spray Gun 테프론 스프레이 건, D.I. Water용

- Recirculating 구조로 정체수에서 발생할 수 있는 박테리아 형성 방지

with OD 1/2" Tubing Set

Cat. No.	Model	Description	EA Price ₩
G09-34-150	421-57	DI Water Gun Kit	

- 구성 : 1/4" FNPT Gun, Recirculating Tee, Male Connector, od1/2" PVC Tubing (2.4 m)
- 재질 : PFA, PTFE
- Entegris

with OD 3/8" Coil Tubing Set

Cat. No.	Model	Description	EA Price ₩
G09-34-250	421-66	DI Water Spray Gun Kit	

- 구성 : PFA Gun, Coiled 1/8" Recirculating Tubing, Flaretek® Recirculating Tee Fitting (3/8" Flaretek x 1/4" NPT x 1/8" Integral Ferrule)
- Entegris

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
G09-34-911	4213211	DI Water Gun Only, 1/4" FNPT	
G09-34-915	21105	Recirculator Tee Assembly	
G09-34-919	603864	Dispense Cap	
G09-34-923	21598	Male Connector	

- Recirculating Tee Assembly
 - For Continuous Flow of DI Water Through the Gun Preventing Possible Bacteria Formation that may Occur in Satagnant Wate
 - od1/8" FEP Tubing 3 m 포함
 - 1/4" MNPT + 1/8" Compression + 1/2" Compression
- Dispense Cap: Ø5.1mm Single Hole, od8.9mm Nozzle Tip
- Male Connetor: 1/4" FNPT Inlet on Gun과 od1/2" Tubing 연결용



Cleanroom Water Gun 클린룸용 워터건, D.I. Water용

Cat. No.	Model	Description	EA Price ₩
G09-695-200	24-0881	Gun Only	201,700
G09-695-205	24-0882	Gun with Blade Hose	248,300
G09-695-210	24-0883	Gun with Blade Hose, Inner Tube & Branch Tee Set	518,000

- 반도체용 세정기에 사용되는 초순수용 워터건
- 액체와 닿는 부위는 모두 PP와 볼스 고무만으로 되어있음
- 울 수지 제품으로 금속 이온의 용출이나 녹 발생의 문제가 없음
- Shower Nozzle 1개 장착 및 Straight Nozzle 1개 포함
- 용도에 따라 본체만, 일반 호스 또는 박테리아 방지 기능 호스를 사용할 수 있음
- 코일 호스는 별도 구매 해야함
- Sansyo

Cat. No.	G09-695-200	G09-695-205	G09-695-210
Model	24-0881	24-0882	24-0883
Temperature	5 ~ 40 °C		
Max. Pressure	0.4 MPa		
Nozzle	Straight, Shower 각 1ea		
Hose	-	Blade Hose	Blade&Inner Hose
Type	본체만	호스 포함	박테리아 방지기능

- Blade Hose : od13.5×h8 mm, Inner Tube : od3×d2 mm

Optional Coil Hose

Cat. No.	Model	길이	재질	EA Price ₩
G09-695-215	24-0884	2.5 m	PP	283,300

- 사용 조건 : 5 - 30°C, Max. 0.2 Mpa
- 1/4" Fitting×2개 포함

