

## Graphite Digestion System / Heating Plate 에서 계속...

Acid Cycle PFA Digestion Vessel Only  
테프론 산순환 포집 분해 용기

Cat. No.	Model	용량	Size (od×h)	EA Price ₩
H04-1995-171	PV-AD-015	15 ml	28.5×100 mm	172,300
H04-1995-174	PV-AD-030	30 ml	28.5×130 mm	179,500
H04-1995-177	PV-AD-060	60 ml	47×145 mm	277,100
H04-1995-180	PV-AD-100	100 ml	47×170 mm	289,200
H04-1995-183	PV-AD-300	300 ml	62×225 mm	330,100

• PD-AD-15 & PV-AD-30 모델은 Spare PFA Cap 포함

## Cap for PV-LP-60 &amp; PV-LP-100



Cat. No.	Model	Description	EA Price ₩
H04-1995-186	PV-CAP-47	PFA Cap Only	48,200

• 60 ml & 100 ml PFA Low Pressure Digestion 용

## Low Pressure PFA Digestion Vessel

Cat. No.	Model	용량	Size (Ø×h)	EA Price ₩
H04-1995-189	PV-LP-015	15 ml	28.5×45 mm	33,700
H04-1995-192	PV-LP-030	30 ml	28.5×75 mm	41,000
H04-1995-195	PV-LP-060	60 ml	47×75 mm	84,300
H04-1995-198	PV-LP-100	100 ml	47×95 mm	95,200

• Bottom - 기본 Flat형 & Round형 변경 가능

High Temp. Heating Plate, SIC Coating  
고온용 흑연 히팅 플레이트, Block 별도

Cat. No.	Model	온도 범위	Plate Size (w×h)	EA Price ₩
H04-1995-101	GDS-HT	Max. 500°C	300×230 mm	5,421,700

- 알루미늄이 아닌 흑연 소재 사용으로 중금속 오염 방지
- SIC 코팅으로 유해 가스에 의한 부식 방지 및 파손 방지
- 블럭형태로 제작하여 분해효율 극대화로 분해시간 단축 (Block 별도 구매)
- 고순도 PFA 테프론 산순환 · 포집 분해용기 사용
- 압력 분해 용기, 일반 비이커등 다양한 용기 호환 가능
- 흑연블럭 사이즈 주문자 방식 생산 가능
- LK Labkorea

## Graphite Block with SIC Coating

Cat. No.	Model	Size (w×d×h, mm)	Type	EA Price ₩
H04-1995-103	GB-SI01	300×230×27, 12홀(Ø56)	Lower block	2,168,700
H04-1995-105	GB-SI02	320×250×18(15), 12홀(Ø56)	Quartz용 Upper block	2,168,700

• 배열 : 3×4열, 12홀

## Quartz Vessel

Cat. No.	Model	용량	Size (Ø×h)	EA Price ₩
H04-1995-110	GB-QV050	50 ml	55×42 mm	201,400

## PTFE Vessel with Cover

Cat. No.	Model	용량	Size (Ø×h)	EA Price ₩
H04-1995-109	GB-PF050	50 ml	53/60×61 mm	397,600

## 상세 사양

Model		GDS-HT
Size	Plate	w300×d230×h25 mm
	Overall	w410×d340×h120 mm
Controller	PID Program Controller	
온도 범위	Max. 500°C (0.1°C 분해능)	
온도 정밀도 / 분포도	±3°C / ±5°C	
Plate Material	SIC Coated Graphite	
Electronic	200 ~ 240V, 50/60hz	

