

Interchangeable heating Block 교환식 히팅 블록, for IKA dry block heater

- IKA Dry Block Heater 전제품에 사용 가능
- 제품 규격은 w95×d76×h51로 동일 (Plate용 제외)
- IKA

for Centrifuge Tubes

Cat. No.	Model	Type	Hole 크기	Hole 수	깊이	Shape	EA Price ₩
H04-05-501	4467600	DB1.1	∅7.9	30	27.6		270,000
H04-05-505	4467700	DB1.2	∅11.5		36.9		270,000
H04-05-509	4467800	DB1.3	∅11.1	20	39.1		270,000
H04-05-513	4467900	DB1.4	∅10.9				190,000
H04-05-517	4468000	DB1.5			38.1		270,000

for Conical Bottom Centrifuge Tubes

Cat. No.	Model	Type	Hole 크기	Hole 수	깊이	Shape	EA Price ₩
H04-05-521	4468100	DB2.1	∅17.1	12	44.5		270,000
H04-05-525	4468200	DB2.2	∅29.0	5	47.6		270,000

for Combination Block

Cat. No.	Model	Type	Hole 크기	Hole 수	Shape	EA Price ₩
H04-05-529	4468300	DB3.1	∅8.3 / 13.9 / 26.2	6/5/3		270,000
H04-05-533	4468400	DB3.2	∅11.1 / 17.1 / 29.0	4/3/2		270,000
H04-05-537	4468500	DB3.3	∅7.9 / 11.1 / 11.5	6/10/5		270,000

• Combination detail

- 4468300 : 6 mm / 12 ~ 13 mm / 25 mm (깊이 : 48.4 / 48.4 / 48.4)
- 4468400 : 1.5 ml / 15 ml / 50 ml (깊이 : 39.1 / 44.5 / 47.6)
- 4468500 : 0.5 ml / 1.5 ml / 2.0 ml (깊이 : 27.6 / 39.1 / 38.1)

for Standard Test Tube

Cat. No.	Model	Type	Hole 크기	Hole 수	깊이	Shape	EA Price ₩
H04-05-541	4468600	DB4.1	∅8.3	30		270,000	
H04-05-545	4468700	DB4.2	∅10.7	24		230,000	
H04-05-549	4468800	DB4.3	∅13.9	16		230,000	
H04-05-553	4468900	DB4.4	∅13.9	20		48.4	250,000
H04-05-557	4469000	DB4.5	∅17.5	12		250,000	
H04-05-561	4469100	DB4.6	∅21.0	8		270,000	
H04-05-565	4469200	DB4.7	∅26.2	6		260,000	
H04-05-569	4469300	DB4.8	∅35.0	4		47.6	280,000
H04-05-573	4469400	DB4.9	∅19.1	12		48.4	280,000

for Standard Test Tube

Cat. No.	Model	Type	Hole 크기	Hole 수	깊이	Shape	EA Price ₩
H04-05-577	4469500	DB5.1	∅12.7	20	30.0		230,000
H04-05-581	4469600	DB5.2	∅15.8		35.0		230,000
H04-05-585	4469700	DB5.3	∅17.8	12	45.0		240,000
H04-05-589	4469800	DB5.4	∅19.7				270,000
H04-05-593	4469900	DB5.5	∅21.7	9	240,000		
H04-05-597	4470000	DB5.6	∅23.8	8	240,000		
H04-05-601	4470100	DB5.7	∅25.8		240,000		
H04-05-605	4470200	DB5.8	∅28.8	6	240,000		
H04-05-609	4470300	DB5.9	∅16.4	15	270,000		

for Cuvettes

Cat. No.	Model	Type	Hole 크기	Hole 수	깊이	Shape	EA Price ₩
H04-05-613	4470700	DB8.1	12.5×6	2	25.4	사진참고	260,000

for Microtitre Plate

Cat. No.	Model	Type	가용 Plate	깊이	EA Price ₩
H04-05-617	4470800	DB7.1	96/384 Well Plate	13.5	330,000

• Dimension of Plate(Double Block)

- w95×d152×h57 : 다른 블록에 2배 크기로 짝수 블록 제품에만 사용 가능

