

Reaction Block for Round Bottom Flask 환저 플라스크용 히팅 교반 블럭

Cat. No.	Model	적용 R-B Flask	Optional Insert	EA Price ₩
H04-132-400	1992-10	100 ml	10, 25, 50 ml	1,184,000
H04-132-410	1992-50	500 ml	200, 250, 300 ml	1,276,900
H04-132-420	1992-60	1 Lit.	500 ml	1,665,400

- Teflon Shield 기본 포함
- Oil Bath나 Heating Mantle 대체용
- 135 mm의 Top Plate를 가진 Hotplate Stirrer용
- Designed to Operate at Temp. up to 180°C
- Chemglass

Optional Insert for 100 ml 용 Reaction Block, H135.4

Cat. No.	Model	적용 R-B Flask	색상	EA Price ₩
H04-132-467	1992-12	10 ml	Gold	237,200
H04-132-475	1992-14	25 ml	Blue	251,400
H04-132-483	1992-16	50 ml	Red	268,600

Optional Insert for 500 ml 용 Reaction Block, H135.5

Cat. No.	Model	적용 R-B Flask	색상	EA Price ₩
H04-132-487	1992-53	200 ml	Turquoise	397,700
H04-132-491	1992-54	250 ml	Clear	397,700
H04-132-495	1992-55	300 ml	Black	397,700

Optional Insert for 1L 용 Reaction Block, H135.6

Cat. No.	Model	적용 R-B Flask	색상	EA Price ₩
H04-132-499	1992-58	500 ml	Purple	606,900



1992-60



사용 예



1992-12



1992-53

Reaction Block for Round Bottom Flask 환저 플라스크용 히팅 교반 블럭, Quarter System

- Chemglass

Block Holder Only, Green Anodized


Cat. No.	Model	Description	EA Price ₩
H04-132-100	1991-P-05	Ø135 mm Hotplate Stirrer 용	237,000

- 135 mm의 Top Plate를 가진 Hotplate Stirrer용
- 4개의 Reaction Block을 거취할 수 있음
- Chemglass

Optional Reaction Block for Round Flask

Cat. No.	Model	적용 R-B Flask	색상	EA Price ₩
H04-132-192	1991-P-30	10 ml	Gold	265,300
H04-132-195	1991-P-32	25 ml	Blue	271,200
H04-132-198	1991-P-34	50 ml	Red	277,300
H04-132-201	1991-P-36	100 ml	Black	283,100

- Ø3 mm Temperature Sensor을 넣기 위한 Hole이 있음



1991-P-36



Hotplate Stirrer는 633 ~ 647 Page 참조

