

## Heating Tape 히팅 테이프

Cat. No.	Model	Width	길이	Watts	EA Price ₩
H04-03-601	MS-HT2506	25.4 mm (1")	0.6 m	200 W	23,000
H04-03-608	MS-HT2510		1 m	380 W	23,000
H04-03-615	MS-HT2515		1.5 m	300 W	34,000
H04-03-622	MS-HT2520		2 m	340 W	45,000
H04-03-629	MS-HT2525		2.5 m	460 W	56,000
H04-03-636	MS-HT2530		3 m	490 W	71,000

- Heating Element는 Glass Fiber로 편직되어 가열, 유연성이 뛰어나
- Max. Temp. : up to 450°C
- Mi-sung



## Stirring Multi Heating Mantle, Analog 교반용 멀티 히팅 맨틀, 3구 or 6구, Flask, EAMDS

- 각기 독립적인 Analog 온도 조절기
- Heating Element는 Glass Fiber로 편직되어 열전달 가열성이 뛰어나
- Ceramic Fiber Insulation 채택으로 사용자 안전을 위한 외부 표면 단열 실현
- Max. Temp. : up to 450°C
- Ø12.7 mm Support Rod 기본 포함
- Timer : 99 hr 59 min
- Stirring Speed : 100 ~ 1200 rpm
- 220 V (요청에 의해 110 V용 공급 가능), 50/60 Hz
- Mi-sung

### 3구 Analog Stirring Multi Heating Mantle

Cat. No.	Model	용량	Watts	EA Price ₩
H04-03-490	MS-EAMDS9601-03	100 ml	80 W×3ea	1,073,000
H04-03-497	MS-EAMDS9602-03	250 ml	150 W×3ea	1,073,000
H04-03-504	MS-EAMDS9603-03	500 ml	260 W×3ea	1,120,000
H04-03-512	MS-EAMDS9704-03	1 L	380 W×3ea	1,190,000

### 6구 Analog Stirring Multi Heating Mantle

Cat. No.	Model	용량	Watts	EA Price ₩
H04-03-518	MS-EAMDS9201-06	100 ml	80 W×6ea	2,122,000
H04-03-525	MS-EAMDS9202-06	250 ml	150 W×6ea	2,122,000
H04-03-532	MS-EAMDS9203-06	500 ml	260 W×6ea	2,210,000
H04-03-539	MS-EAMDS9204-06	1 L	380 W×6ea	2,319,000



## Digital Temperature Indicator 디지털 온도 표시기

Cat. No.	Model	Temp. Range	EA Price ₩
H04-03-811	MS-HTI610K	K-type (0 ~ 1200°C)	198,000
H04-03-812	MS-HTI610Pt	PT-type (-99.9 ~ 199.9°C)	198,000

- 각 포인트 온도 표시 자동 / 수동 전환 가능
- Point : 5개
- 자동 시간 조절기능 (1 ~ 10 sec)
- Size : w205×d220×h70 mm
- AC 220V, 50/60 Hz
- Sensor는 별도 구매
- Mi-sung

