

Biopsy Embedding Cassette 바이옵시 임베딩 카세트

Cat. No.	Model	색상	PK Price ₩		CS Price ₩	
H06-660-171	0109-1100-16	White	250ea	25,100	8pk	187,100
H06-660-175	0109-1100-17	Yellow	250ea	25,100	8pk	187,100


- 재질 : POM, 100% Resistant to Histological Solvents
- Square Hole : w0.9×d0.9 mm
- 2 Side Labelling Area, Front 45 ° Angle
- LK Labkorea



Tissue Embedding Cassette 티슈 임베딩 카세트, One-piece

Cat. No.	Model	색상	PK Price ₩		CS Price ₩	
H06-660-179	0106-1100-16	White	250ea	24,600	8pk	183,700
H06-660-183	0106-1100-17	Yellow	250ea	24,600	8pk	183,700


- 재질 : POM, 100% Resistant to Histological Solvents
- One-piece & Snap-latch Design
- Rectangle Hole : w5.0×d1.0 mm
- 2 Side Labelling Area, Front 45 ° Angle
- LK Labkorea



Tissue Embedding Cassette 티슈 임베딩 카세트, Two-pieces

Cat. No.	Model	색상	PK Price ₩		CS Price ₩	
H06-660-187	0100-1101-16	White	500ea	52,900	4pk	201,700
H06-660-191	0100-1101-17	Yellow	500ea	52,900	4pk	201,700

- 재질 : POM, 100% Resistant to Histological Solvents
- Two-pieces Design, Cassette와 Lid 분리형
- Rectangle Hole : 5.0 mm
- 2 Side Labelling Area, Front 45 ° Angle
- LK Labkorea



Biopsy Foam Pad 임베딩 패드

Cat. No.	Model	Size	색상	PK Price ₩
H06-660-195	0110-0001	w25×d32×h2 mm	White	1000ea 29,500

- 재질 : PU Foam
- Biopsy를 제자리에 고정시키는 제품
- LK Labkorea



Stainless Steel Base Mold 스테인레스 베이스 몰드

Cat. No.	Model	Net Well Size	PK Price ₩
H06-660-199	0203-0001	w7×d7×h7 mm	10ea 13,000
H06-660-203	0203-0002	w15×d15×h7 mm	10ea 13,000
H06-660-207	0203-0003	w24×d24×h7 mm	10ea 13,000
H06-660-211	0203-0004	w37×d24×h7 mm	10ea 13,000
H06-660-215	0203-0005	w30×d24×h7 mm	10ea 13,000

- 재질 : Stainless Steel 304
- Size : w35×d52×h11 mm
- Fine Polished Surface & Well-Grounded Corner
- LK Labkorea

