

Disposable Base Mold 일회용 베이스 몰드

Cat. No.	Model	Net Well Size	PK Price ₩	CS Price ₩
H06-660-219	0203-0006	w7×d7×h5 mm	500ea 34,000	10pk 317,200
H06-660-223	0203-0007	w17×d17×h5 mm	500ea 34,000	10pk 317,200
H06-660-227	0203-0008	w25×d25×h5 mm	500ea 34,000	10pk 317,200
H06-660-231	0203-0009	w25×d32×h5 mm	500ea 34,000	10pk 317,200
H06-660-235	0203-0010	w35×d40×h5 mm	500ea 34,000	10pk 317,200

- 재질 : PVC
- Size : w37×d50×h12 mm
- LK Labkorea



Block Storage Cabinet 블럭 저장 박스

Cat. No.	Model	서랍 수	보관 Block 수	EA Price ₩
H06-660-239	0510-0040	6개	up to 21000개	201,300


- 재질 : Heavy-duty ABS
- Size : w245×d445×h400 mm
- 여러개의 Cabinet을 중첩하여 쌓아 올릴 수 있음
- Plastic Covered Label 포함
- LK Labkorea



Specimen Storage Cabinet 시료 보관 캐비닛 모듈

Cat. No.	Model	Module	높이	EA Price ₩
H06-660-243	0510-0006	Base	76 mm	43,600
H06-660-247	0510-0007	Lid	26 mm	35,300
H06-660-251	0510-0008	Slide Storage	126 mm	166,000
H06-660-255	0510-0009	Paraffin Block Storage	91 mm	151,500

- 재질 : White Coated Steel
- Size : w485×d480 mm
- Specimen Filling & Storage용
- Base와 Lid 사이에 용도에 맞추어 개별 Storage Rack을 추가하는 모듈형
- Slide Storage (0510-0008) : 14 Individual Drawers
(각 최대 5600 Slide 보관 가능)
- Paraffin Block Storage (0510-0009) : 14 Individual Drawers
(각 최대 900 Block 보관 가능)
- LK Labkorea



Paraffin Wax, HistoPlast® 파라핀 왁스

Cat. No.	Model	Melting Point	EA Price ₩
H06-660-167	0200-0007	56 ~ 58°C	11,800

- Compound of Purified Paraffin and Plastic Polymers
- without Dimethyl Sulfoxide (DMSO)
- Delivers Smooth Cutting, Minimized Compression & Cell Disruption
- Cutting Thickness : low to 2 um
- Convenient Pearl Shape
- LK Labkorea

