

IKA Homogenizer, ULTRA-TURRAX® 소용량 균질기, T10 basic

Cat. No.	Model	Type	작업 가능 용량	EA Price ₩
H01-05-205	3737000	T10 basic	0.5 ~ 100 ml	3,140,000
H01-05-125	3737000P	T10 Package		4,357,500
H01-05-920	3372000	R200	Clamp	130,000
H01-05-935	2437700	H44	Boss Head Clamp	60,000
H01-05-960	3386000	R104	Stand	280,000

- 본체만의 가격이며, Tool 및 기타 Accessory는 별도
- T10 Package : T10 basic, S10N-8G, H44, R104, R200
- Hand-held 및 Stand를 이용하여 사용 가능
- Dispersing elements를 쉽게 교체할 수 있는 커플링(Coupling)방식
- IKA

Type	T10 basic	
Motor Rating in/out	125W / 75W	
Volume Range	1 ~ 100 ml (H2O)	
Max. Viscosity	5,000 mPas	
Speed Range	8,000 ~ 30,000 rpm 8,000 ~ 26,000 rpm	
Speed 조절	6 단 (Scale)	
Size	Extension Arm	∅8×L100 mm
	Overall	w63×d57×h210 mm
Power	230 V, 50/60 Hz	

Tool / Dispersing Element for T10 basic

Cat. No.	H01-05-514	H01-05-520	H01-05-526
Model	3304000	3305500	3370100
Type	S10N-5G	S10N-8G	S10N-10G
Working Range	0.5 ~ 10 ml	1 ~ 50 ml	1 ~ 100 ml
Stator Diameter	5 mm	8 mm	10 mm
Rotor Diameter	3.8 mm	6.1 mm	7.6 mm
Immersion min/max	20 / 75 mm	20 / 95 mm	20 / 100 mm
Shaft length	92 mm	115 mm	
Ultimate Suspension	5 - 25 um		
Ultimate Emulsion	1 - 10 um		
EA Price ₩	710,000	960,000	1,010,000

Plastic Dispersing Element for T10 basic

Cat. No.	H01-05-817	H01-05-829
Model	3433225	3433225
Type	S10D-7G-KS-65	S10D-7G-KS-110
Working Range	1 ~ 20 ml	1 ~ 40 ml
Stator Diameter	7 mm	
Rotor Diameter	4.8 mm	
Immersion min/max	20 / 50 mm	20 / 90 mm
Shaft length	65 mm	110 mm
Max. Temperature	100°C, Autoclavable	
재질	Polycarbonate (PC) & Polysulfon (PCU), FDA 승인	
PK Price ₩	24ea 520,000	24ea 580,000

- Ideal for Homogenizing Tissue Sample

Sawtooth (ST) dispersing Tool for T10 basic

Cat. No.	H01-05-831	H01-05-833
Model	4446500	4446700
Type	S10N-8G-ST	S10N-10G-ST
Working Range	1 ~ 50 ml	1 ~ 100 ml
Stator Diameter	8 mm	10 mm
Rotor Diameter	6.1 mm	7.5 mm
Immersion min/max	22 / 95 mm	23 / 100 mm
Shaft length	115 mm	
Rotor/Stator Gap	0.25 mm	0.25 mm
재질	PTFE, AISI316L	
EA Price ₩	1,240,000	1,310,000

