


## IKA Homogenizer, ULTRA-TURRAX® 계속...

## Tool/Dispersing Element for T25 digital



적용 Model	T25 digital				
Cat. No.	H01-05-544	H01-05-550	H01-05-556	H01-05-562	H01-05-568
Model	1024200	0594000	0593400	1713300	1713800
Type	S25N-8G	S25N-10G	S25N-18G	S25N-25G	S25N-25F
Working Range	1 ~ 50 ml	1 ~ 100 ml	10 ~ 1500 ml	50 ~ 2000 ml	100 ~ 2000 ml
Stator Diameter	8 mm	10 mm	18 mm	25 mm	25 mm
Rotor Diameter	6.1 mm	7.5 mm	12.7 mm	17 mm	18 mm
Immersion min/max	27 / 85 mm	22 / 85 mm	40 / 165 mm	40 / 165 mm	40 / 165 mm
Shaft Length	108 mm	105 mm	194 mm	194 mm	194 mm
Ultimate Suspension	10 - 50 um			15 - 50 um	5 - 25 um
Ultimate Emulsion	1 - 10 um				1 - 5 um
Max. Temperature	180 °C				
재질	PTFE, AISI 316L				
EA Price ₩	1,420,000	1,440,000	1,870,000	1,780,000	1,750,000



적용 Model	T25 digital				
Cat. No.	H01-05-855	H01-05-858	H01-05-861	H01-05-865	H01-05-877
Model	2348000	2466900	2404000	3452200	3452100
Type	S25KV-18G	S25KV-25G	S25KV-25F	S25D-10G-KS	S25D-14G-KS
Working Range	10 ~ 1500ml	50 ~ 2000ml	100 ~ 2000ml	10 ~ 100ml	10 ~ 500ml
Sator Diameter	18 mm	25 mm	25 mm	10 mm	14 mm
Rotor Diameter	12.7 mm	17 mm	18 mm	6.75 mm	9.5 mm
Immersion min/max	40 / 225 mm	40 / 225 mm	40 / 225 mm	15 / 85 mm	15 / 85 mm
Shaft Length	270 mm	270 mm	270 mm	150 mm	150 mm
Ultimate Suspension	10 - 50 um	15 - 50 um	5 - 25 um	10 - 50 um	10 - 50 um
Ultimate Emulsion	1 - 10 um	1 - 10 um	1 - 5 um	5 - 20 um	5 - 20 um
Max. Temperature	220 °C			100 °C	
재질	FFPM / SIC, AISI 316L			Polycarbonate (PC) & Polyetheretherketon (PEEK)	
EA Price ₩	2,770,000	2,830,000	2,970,000	PK Price ₩ 12ea 690,000	PK Price ₩ 12ea 690,000



적용 Model	T25 digital			
Cat. No.	H01-05-879	H01-05-881	H01-05-882	H01-05-884
Model	4446900	4447100	4447300	4447500
Type	S25N-8G-ST	S25N-10G-ST	S25N-18G-ST	S25N-25G-ST
Working Range	1 ~ 50 ml	1 ~ 100 ml	10 ~ 1500 ml	50 ~ 2000 ml
Stator Diameter	8 mm	10 mm	18 mm	25 mm
Rotor Diameter	6.1 mm	7.5 mm	12.7 mm	17 mm
Immersion min/max	27 / 85 mm	22 / 85 mm	40 / 165 mm	40 / 165 mm
Shaft length	108 mm	105 mm	194 mm	194 mm
재질	PTFE / AISI 316L			
EA Price ₩	1,910,000	1,920,000	2,010,000	2,030,000

• 그 외 S25KD-Series는 홈페이지 참조