

## IKA Homogenizer for Large Volume 대용량 균질기, T50 digital

Cat. No.	Model	Type	용량	EA Price ₩
H01-05-220	3787000	T50 digital	0.25 ~ 30 L	2,508,000

- 본체만의 가격이며, Tool 및 기타 Accessory는 별도
- 속도표시 : LED Digital Display
- 230V, 50/60 Hz
- IKA

Cat. No.	H01-05-220
Type	T50 digital
Motor Rating in/out	1,100W / 700W
Volume Range	0.25 ~ 30 L
Max. Viscosity	5,000 mPas
Speed Range	600 ~ 10,000 rpm
Extension Arm	Ø16×L220 mm
Size	w115×d139×h355 mm

### Optional Clamp

Cat. No.	Model	Type	Stand Rod	Extension Arm	EA Price ₩
H01-05-945	2657800	R270	Ø25-36 mm	Ø5 ~ 21 mm	260,000
H01-05-950	2664000	R271	Ø34 mm	Ø16 mm	260,000

### Optional Stand

Cat. No.	Model	Type	높이 조절	Size	EA Price ₩
H01-05-985	1412100	R2723	O	34×620~1010 mm	1,800,000



T50 digital

### Tool/Dispersing Element for T50 digital

Cat. No.	H01-05-574	H01-05-580	H01-05-586
Model	8003000	8003300	8003900
Type	S50N-G45G	S50N-G45M	S50N-G45F
Working Range	0.5 ~ 20 L	0.5 ~ 15 L	0.25 ~ 10 L
Stator Diameter	45 mm		
Rotor Diameter	36 mm	40.5 mm	40 mm
Immersion min/max	70 / 250 mm		
Shaft length	300 mm	290 mm	
재질	AISI 316L, PTFE bearing (내열 온도 180°C)		
Ultimate Suspension	40 ~ 100 um	25 ~ 50 um	10 ~ 30 um
Ultimate Emulsion	10 ~ 30 um	5 ~ 20 um	1 ~ 10 um
EA Price ₩	3,440,000	3,440,000	3,430,000



S50N-G45G



S50N-G45M



S50N-G45F

### Special Dispersing Element for T50 digital



Cat. No.	H01-05-608	H01-05-614	H01-05-620	H01-05-626	H01-05-632
Model	1689300	8006300	8005100	1289800	1243300
Type	R50	S50N-W80SMK	S50N-W65SK	R1405	R1402
Working Range	0.25 ~ 30 L	1 ~ 50 L	1 ~ 10 L	0.25 ~ 30 L	1 ~ 30 L
Rotor Diameter	50 mm (R1402포함)	80 mm	65 mm	45 mm	42 mm
규격	Stirring Shaft	Jet Mixer Head	Cutting Head	"R50"용 dissolver	
Immersion	180 mm	140 / 350 mm	80 / 350 mm	Optional	"R50"에 기본 포함
재질	AISI 316L				
EA Price ₩	4,200,000	3,610,000	12,050,000	440,000	580,000

- 그 외 특수한 Tool은 홈페이지 참조