

I Weighing Bottle 평량병, LukeGL®

- Borosilicate Glass 3.3
- LukeGL®

Short Form

Cat. No.	Model	용량	Size (Ø×h)	Joint	EA Price ₩
B08-96-010	70-1010	7 ml	30×30 mm	24/12	18,400
B08-96-020	70-1020	18 ml	35×40 mm	29/12	19,000
B08-96-030	70-1030	23 ml	40×40 mm	34/12	20,000
B08-96-040	70-1050	54 ml	50×50 mm	45/12	20,400
B08-96-050	70-1100	102 ml	60×60 mm	50/12	25,200

Tall Form

Cat. No.	Model	용량	Size (Ø×h)	Joint	EA Price ₩
B08-96-070	71-1010	9 ml	25×40 mm	19/12	18,200
B08-96-080	71-1020	19 ml	30×50 mm	24/12	18,700
B08-96-090	71-1030	26 ml	35×50 mm	29/12	19,300
B08-96-100	71-1050	45 ml	40×60 mm	34/12	20,400
B08-96-110	71-1100	105 ml	50×80 mm	45/12	21,200



Short Form

Tall Form

I Tall Cylindrical Weighing Bottle, Kimble 장형 평량병

Cat. No.	Model	용량	EA Price ₩	PK Price ₩
B17-137-001	15146-2540	12 ml	50,800	6ea 284,900
B17-137-003	15146-2550	16 ml	51,300	6ea 287,600
B17-137-005	15146-3060	30 ml	57,900	6ea 324,600
B17-137-007	15146-4050	45 ml	70,700	6ea 396,400
B17-137-009	15146-4080	70 ml	73,700	6ea 413,000
B17-137-011	15146-40100	92 ml	82,500	6ea 462,400

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- DWK Lifescience / Kimble®

상세규격

Model	Body (ID×h)	전체 높이	Joint	Approx. Weight
15146-2540	25×40 mm	52 mm	29/12	20 g
15146-2550	25×50 mm	62 mm		23 g
15146-3060	30×60 mm	72 mm	34/12	32 g
15146-4050	40×50 mm	64 mm		55 g
15146-4080	40×80 mm	94 mm	45/12	72 g
15146-40100	40×100 mm	114 mm		83 g



※ 15146-4080은 1×3" Micro TLC Plate와 같이
TLC Chamber로 사용하기에 알맞음

I Aluminum Weighing Bottle 알루미늄 평량병

Cat. No.	Model	용량	Size	EA Price ₩
B19-125-005	AWB020	20 ml	Ø50×h20 mm	5,500
B19-125-010	AWB025	25 ml	Ø60×h20 mm	6,300
B19-125-015	AWB075	70 ml	Ø60×h40 mm	7,000

- 표면처리가 된 소용량의 평량용
- LK Labkorea



I Plastic Weighing Bottle 플라스틱 평량병

Cat. No.	Model	용량	Size	EA Price ₩
B17-339-111	PWB-020	20 ml	Ø32×h40 mm	2,050
B17-339-115	PWB-030	30 ml	Ø32×h65 mm	2,170
B17-339-117	PWB-050	50 ml	Ø40×h52 mm	2,530
B17-339-121	PWB-070	70 ml	Ø52×h65 mm	2,890

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- LabPlasti®

