

Ductless ArmHood 이동형 필터식 암후드

Cat. No.	Model	포함 Filter	용도	EA Price ₩
H02-260-201	ULT_JUMBO	Carbon	가스, 약취, 증기	4,891,500
H02-260-205	ULT JUMBO T2	HEPA	분진, 스모크	5,301,200
H02-260-209	ULT JUMBO T3	Carbon + HEPA	가스, 약취, 증기 분진, 스모크	5,072,200

- 사용가능 환경 : 납땜 및 레이저 작업시 연기 제거용, 증기와 가스 제거용, 크고 작은 분진 제거용
- 사용 환경에 따라 필터 변경 가능
- 컴팩트한 사이즈로 가볍고 저소음이며 이동이 편리함
- 흡입구 판과 흡입구 관의 위치 변경 가능 (최적화된 작업 환경)
- 배기 노즐이나 호스를 사용하여 배출량 조절이 가능함
- Ø50mm, Ø63mm, Ø75mm 흡입구 사용 가능
- ULT JUMBO용 Arm 포함
- 디지털 컨트롤이 가능하며 다른 기기와 쉽게 연결할 수 있음 (Option)
- 핸들 부착형 리모트 컨트롤 (Option)
- LK Labkorea

Filter for ULT JUMBO

Cat. No.	Model	Filter Type	용도	EA Price ₩
H02-260-570	JUMBO FC	Carbon	가스, 약취, 증기 제거용	518,000
H02-260-572	JUMBO FH	HEPA	분진, 스모크 제거용	783,100
H02-260-576	JUMBO FCH	Carbon+HEPA	가스, 약취, 증기, 분진, 스모크 제거용	698,700
H02-260-580	JUMBO FA	100% Chemisorption (산 전용)	산 제거용	1,277,100
H02-260-584	JUMBO FCA	50% Carbon + 50% Chemisorption	가스, 약취, 증기, 산 제거용	927,700



Model	ULT_JUMBO
Max. Flow Rate	170 m³/h
Max. Vacuum	2800 Pa
Power Output (flow rate)	80 m³/h at 1400 Pa
Electrical Connection	230V / 50 Hz
Power Input	150 W
Size	w450 × d340 × h530 mm
Weight	18 kg
Noise Level	49 dB(A)
Material	Impact-resistant ABS resistant to chemicals resistant to corrosion adjustable blow-out direction optional air outlet
Arm reach	1125 mm / 3 joint / Ø50 mm

필터식 시약장은 528 ~ 529 Page 참조

