

Jumbo Hotplate 대형 핫 플레이트

PID Control

Cat. No.	Model	Plate Size	EA Price ₩
H03-04-050	DLP250	250×250 mm	888,900
H03-04-060	DLP300	300×300 mm	908,100
H03-04-080	DLP430	400×300 mm	1,033,600
H03-04-090	DLP630	600×300 mm	1,147,300

Program Control

Cat. No.	Model	Plate Size	EA Price ₩
H03-04-055	DLP250P	250×250 mm	1,123,700
H03-04-065	DLP300P	300×300 mm	1,147,300
H03-04-087	DLP430P	400×300 mm	1,305,300
H03-04-095	DLP630P	600×300 mm	1,450,100

- 일반 가열, 적정(Titration), Wafer Baking용으로 적합
- 우수한 온도 분포도와 Plate 평탄도 • Plate 최대 하중무게 : 5kg
- Anodizing 가열판을 사용하여 내열성, 내약품성 우수
- Graphic LCD Display, 간편하고 직관적인 Controller 조작
- Auto-tuning, Temperature Calibration 기능
- 2가지 타입의 타이머. 구동과 동시에 시작 / 목표 온도 도달 후 시작
- Plate 주변 보호 마스크가 있어 사용자가 안전하고 화재 위험이 적음
- Automatic Fault Sensing System 기기 스스로 이상 상황을 감시 (특히 등류)
- 사용자의 요청에 의해 Program Control Model 공급 가능
- CE 인증 모델은 요청에 의하여 공급 가능하며 Over Temp. Limiter 추가
- Program Control
 - 온도 Overshoot 억제 및 정밀 제어 가능
 - 2 Pattern, 18 Segment 프로그램 설정 가능
- Validation Service (IQ, OQ, PQ)
- PL Insurance (제조물 배상 책임 보험)
- LK Labkorea

Validation Service

Cat. No.	Model	Description	EA Price ₩
H03-91-999	VA-JH100	Validation for Jumbo Hotplate	

- 당사 기본 양식에 따른 IQ/OQ Validation Service 비용이며, 지역이나 고객 요청사항에 따라 비용이 변경될 수 있음

Cat. No.	H03-04-050 / H03-04-055	H03-04-060 / H03-04-065	H03-04-080 / H03-04-087	H03-04-090 / H03-04-095
Model	DLP250 / DLP250P	DLP300 / DLP300P	DLP430 / DLP430P	DLP630 / DLP630P
Plate Size (w×d)	250×250 mm	300×300 mm	400×300 mm	600×300 mm
Controller	Control			
	PID Control, Auto-tuning / Program Control			
	Display			
	GLCD (Graphic LCD)			
Temperature	Resolution			
	0.1℃			
	Timer			
	Run start / Wait start (99hr 59min)			
Temperature	Range			
	주변온도 +10℃ to +320℃			
	Accuracy at 100℃			
±0.8℃				
Temperature	Uniformity at 100℃			
	±1.5℃	±1.7℃	±1.9℃	±2.4℃
Size (w×d×h, mm)	258×258×150	308×308×150	408×308×150	608×308×150
Weight	6.7 kg	9.5 kg	12.5 kg	15.3 kg
Material	Plate			
	Anodized Aluminium			
Material	Exterior			
	Powder Coated Steel			
Electric Supply	Power			
	1 Phase, 220 Vac, 50/60 Hz			
	Max Consumption			
1.2 kW (5.5 A)	1.5 kW (6.8 A)	1.8 kW (8.2 A)	2.4 kW (11 A)	
Power Line				
Standard Plug				

Jumbo Hotplate, "CE" 대형 핫 플레이트, CE 인증

PID Control

Cat. No.	Model	Plate Size	EA Price ₩
H03-04-010	DLP250CE	250×250 mm	1,154,500
H03-04-020	DLP300CE	300×300 mm	1,188,100
H03-04-030	DLP430CE	400×300 mm	1,305,900
H03-04-040	DLP630CE	600×300 mm	1,416,800

- Over Temp. Limiter 추가 외의 모든 내용은 상단과 동일

Program Control

Cat. No.	Model	Plate Size	EA Price ₩
H03-04-015	DLP250PCE	250×250 mm	1,405,100
H03-04-025	DLP300PCE	300×300 mm	1,446,000
H03-04-035	DLP430PCE	400×300 mm	1,588,100
H03-04-045	DLP630PCE	600×300 mm	1,724,700

- LK Labkorea

