

Digital High Temperature Hotplate, Max. 550°C 고온용 디지털 핫플레이트, 184×184 Plate

Cat. No.	Model	Max. Temp	Plate Size	EA Price ₩
H03-238-010	HP550-S	550°C	w184×d184 mm	580,800
H03-238-015	HP550-S-PK2		HP550-S Package	540,200

- LED Digital Display : 설정온도 디지털 표시
- Glass Ceramic Plate로 내화학성 및 우수한 열 전도성
- 최대 550°C까지 견딜 수 있는 안정성
- Safety Circuit System : Cut off Temperature 580°C
- Off 상태에서 50°C 초과 시에 "HOT" 경고 알람
- 외부 온도 센서(PT1000)를 통해 정확한 온도제어 가능
- Protection class acc. to DIN EN 60529 (IP21)
- DLAB

Cat. No.		H03-238-010
Model		HP550-S
Plate Size (w×d)		184×184 mm
Temperature	Range	Max. 550°C
	Increment	5°C
	Display Accuracy	±1°C
Controller Accuracy		±0.2°C(PT1000), ±10°C(Working Plate)
Heating Output / Power		1000 W / 1010 W
Heating Warning		50°C
Plate Material		Glass Ceramic
Size		w215×d360×h112 mm
Power		200-240V, 50/60 Hz

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
H03-238-572	D-111000	Temp. Sensor (PT1000)	105,500
H03-238-575	D-111001	Support Kit	68,000

* Support Kit 구성품 : Support Rod, Boss Head Clamp, Holding Rod



- 소형 핫플레이트는 637 Page 참조



Large High Temperature Hotplate, Max. 500°C 대형 고온용 핫플레이트, 254×254 Plate

Cat. No.	Model	Max. Temp	Plate Size	EA Price ₩
H03-238-020	HP500-Pro	500°C	w254×d254 mm	797,000

- LCD Digital Display : 설정온도 및 작동시간 디지털 표시
- 10인치의 Glass Ceramic Plate로 내화학성 및 우수한 열 전도성
- 최대 500°C까지 견딜 수 있는 안정성
- Safety Circuit System : Cut off Temperature 540°C
- Off 상태에서 50°C 초과 시에 "HOT" 경고 알람
- 외부 온도 센서(PT1000)를 통해 정확한 온도제어 가능
- SD Card를 이용하여 작동중일 때도 데이터를 실시간으로 기록할 수 있음
- Protection class acc. to DIN EN 60529 (IP21)
- DLAB

Cat. No.		H03-238-020
Model		HP550-Pro
Temperature	Display	LCD
	Range	Max. 550°C
	Increment	1°C 단위
	Display Accuracy	±0.1°C
Controller Accuracy		±0.2°C(PT1000), ±1°C(Working Plate <100°C), ±1%(Working Plate >100°C)
Heating Output / Power		1350 W / 1400 W
Heating Warning		50°C
Plate Material		Glass Ceramic
Size		w328×d424×h138 mm
Power		200-240V, 50/60 Hz

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
H03-238-572	D-111000	Temperature Sensor (PT1000)	105,500
H03-238-575	D-111001	Support Kit	68,000

