

Digital Hotplate Stirrer with Stainless Steel Plate 디지털 가열 자력 교반기, MS-H-Pro^T or MS-H-Pro⁺

Cat. No.	Model	용량	Plate Size	Timer	EA Price ₩
H03-238-217	MS-H-Pro ^T	20 L	Ø135 mm	Yes	640,500
H03-238-219	MS-H-Pro ^T Pk2	MS-H-Pro ^T Package			697,300
H03-238-211	MS-H-Pro ⁺	20 L	Ø135 mm	Yes	583,500
H03-238-215	MS-H-Pro ⁺ Pk2	MS-H-Pro ⁺ Package			640,500

- Digital Type의 가열 자력 교반기
- Plate 재질 : Ceramic Coated Stainless Steel
- MS-H-Pro^T 모델은 Timer 기능이 내장되어 있음
- Temp. Sensor를 연결하여 외부시료의 정확한 온도제어가 가능
- Bushing according to DIN 12878 for the use of contact thermometers
- Package 구성 : 본체, Temp. Sensor, Support Kit
- DLAB

Model	MS-H-Pro ^T	MS-H-Pro ⁺
Timer Function	Yes (1 min to 99 hr 59 min)	No
Stirring Quantity (H ₂ O)	20 L	
Motor Rating	Input 18 W / Output 10 W	
Speed Range / Resolution	100 ~ 1500 rpm / ±1 rpm	
Heating Output / Rate (H ₂ O)	500 W / 6 k/min	
Temp. Range	RT~340°C	
Temp. Accuracy	±1 °C (Control Accu. w/ Sensor : ±0.2 °C)	
Display	LCD (Temperature & Speed)	
Adjustable Safety Circuit	100 ~ 350°C	
Plate Size	Ø135 mm (Ceramic Coated Stainless Steel 재질)	
Data Connection	RS232	
"Hot" Warning	50°C	
Size	w160×d280×h85 mm (2.8 kg)	
Power	230V, 50/60 Hz, 550 W	

Spare Accessory

Cat. No.	Model	Description	EA Price ₩
H03-238-572	D-111000	Temperature Sensor (PT1000)	105,500
H03-238-575	D-111001	Support Kit	68,000

- Support Kit 구성품 : Support Rod, Boss Head Clamp, Holding Rod



Body only

CE



MS-H-Pro^T / MS-H-Pro⁺

Digital Hotplate Stirrer with Aluminum Plate 디지털 가열 자력 교반기, MS-H-Pro^A

Cat. No.	Model	용량	Plate Size	Timer	EA Price ₩
H03-238-230	MS-H-Pro ^A	20 L	Ø135 mm	No	626,700
H03-238-232	MS-H-Pro ^A Pk2	MS-H-Pro ^A Package			792,000

- Plate 재질 : Ceramic Coated Aluminum Plate (높은 열전도율)
- Temp. Sensor를 연결하여 외부시료의 정확한 온도제어가 가능
- Bushing according to DIN 12878 for the use of contact thermometers
- 3 Heating Mode : Fast Heating, Standard Heating, Stable Heating
- Heating Power : 600W
- Ø135 mm 원형 Heating Block 과 사용에 알맞음
- 기타 자세한 사항은 상부 MS-H-Pro⁺ 와 동일 (단, Timer Function 제외)
- Package 구성 : 본체, Temp. Sensor, Support Kit
- DLAB

Spare Accessory

Cat. No.	Model	Description	EA Price ₩
H03-238-572	D-111000	Temperature Sensor (PT1000)	105,500
H03-238-575	D-111001	Support Kit	68,000

- Support Kit 구성품 : Support Rod, Boss Head Clamp, Holding Rod



MS-H-Pro^A

