

## Large Hotplate Stirrer, 254 × 254 Plate 대형 디지털 고온용 가열 자력 교반기, up to 500°C

Cat. No.	Model	용량	Plate Size	EA Price ₩
H03-238-700	MS10-H500-Pro	30 L	254 x 254 mm	996,100
H03-238-715	MS10-H500-Pro PK 2	MS10-H500-Pro Package		1,038,900

- 고온용 Glass Ceramic Plate 장착
- 10인치 크기의 최대 30L 까지 수용 가능한 대형 가열 자력 교반기
- 외부 센서 (PT1000)를 연결하여 시료의 온도 조절 가능 (±0.2°C)
- Hot Indicator Warning Lamp 장착 ( 50°C 이상시 )
- SD Card를 이용하여 작동중일 때 데이터를 실시간으로 기록할 수 있음
- Package 구성 : 본체, Temp. Sensor, Support Kit
- DLAB

### Spare Accessory

Cat. No.	Model	Description	EA Price ₩
H03-238-572	D-111000	Temperature Sensor (PT1000)	105,500
H03-238-575	D-111001	Support Kit	68,000

- Support Kit 구성품 : Support Rod, Boss Head Clamp, Holding Rod

Model	MS10-H500-Pro
Stirring Quantity(H <sub>2</sub> O)	30 L
Plate Size	254 x 254 mm ( 10 inch )
Display	LCD
Speed / Resolution	100 ~ 1500 rpm,
Temp. Range	RT-500°C (1°C 단위 조절)
Temp. Control Accuracy	±1°C (<100°C), ±1% (>100°C)
Temp. Display Accuracy	±0.1°C
External Temp. Sensor	PT1000 (Accuracy ±0.2°C)
Heating Output	1350 W
"Hot" Warning	50°C
Data Connection	SD card
Protection Class	IP21
Size, w x d x h	328 x 424 x 138 mm / 7.2 kg
Power	200-240v / 50/60 Hz / 1400w



CE

Package 구성 모습



CE

MS10-H500-Pro

## Digital Hotplate Stirrer 소형 디지털 가열 자력 교반기, up to 380°C

Cat. No.	Model	용량	Type	EA Price ₩
H03-238-611	MS-H380-Pro	5 L	Hotplate Stirrer	355,800
H03-238-615	HP380-Pro		Hotplate	298,900
H03-238-617	MS6-Pro		Magnetic Stirrer	284,700
H03-238-613	MS-H380-Pro Pk2	MS-H380 Pro Package		547,500

- Plate : 140 x 140 mm, Ceramic Coated Aluminum
- Brushless DC Motor를 사용 정확한 Speed Control
- 외부 센서를 연결하여 시료의 온도 조절 가능 (±0.5°C)
- Safety Circuits : Cut Off Temperature 420 °C
- Protection class acc. to DIN EN 60529 (IP21)
- Package 구성 : 본체, Temp. Sensor, Support Kit
- DLAB

Model	MS-H380-Pro
Stirring Quantity(H <sub>2</sub> O)	5 L
Plate Size	140 x 140 mm ( Ceramic Coated Aluminum )
Display	LCD, Temperature & Speed
Speed	200 ~ 1500 rpm,
Temp. Range	RT-380°C, 1°C 단위 조절
Temp. Display Accuracy	±0.1°C
External Temp. Sensor	PT1000 (Accuracy ±0.5°C)
Size, w x d x h	320x180x108 mm / 2.2 kg
Power	510 W, 200-240 V, 50/60 Hz



CE

MS-H380-Pro



MS6-Pro