

## Accessory for IKA Hotplate Stirrer / Hotplate IKA 가열 자력 교반기 / 가열판 악세서리

### Temperature Controller

Cat. No.	Model	Type	Range	pH	EA Price ₩
H03-05-506	3378000	ETS-D5	-50 ~ 450°C	-	690,000
H03-05-512	3378100	ETS-D6		0 ~ 14	660,000

- PID Control과 Reset 기능
- Stainless Steel Sensor H62 기본 포함
- PH 값 측정 모듈 내장형 (ETS-D6 : Electrode는 별도임, BNC Type)
- Protection class acc. to DIN EN 60529 (IP54)
- HS Series, RCT basic, RET basic, RH basic과 연결 가능
- IKA

Cat. No.	H03-05-506	H03-05-512
Type	ETS-D5	ETS-D6
Temperature	Range	-50 ~ 450°C
	Resolution	0.1 K
	Accuracy	0.1 K
	Deviation	±0.5 K
pH	Range	0 ~ 14 pH
	Resolution	±0.01 pH
	Accuracy	±0.1 pH
	pH Connection	BNC Bushing
Plug	6 pin DIN 45322	
Connection	DIN 12878 Class 2	
Supply Voltage	8 ~ 16 VDC	
Size	w82 × d22 × h83 mm	w96 × d45 × h98 mm

### Spare Temperature Sensor for "ETS-D" Series

Cat. No.	Model	Type	Description	EA Price ₩
H03-05-518	2735451	H62.51	Stainless Steel	240,000
H03-05-530	2735600	H70	Extension Cable, 1 m	440,000

- Size : Ø3 × L230 mm

### RH Basic & RH Digital용 온도 센서

Cat. No.	Model	Type	재질	EA Price ₩
H03-05-524	2735551	H66.51	Stainless Steel, Glass-coated	630,000
S09-05-500	4499900	H66.53	Stainless Steel, SafeCoat	500,000

- 2735551 Size : Ø6 × L260 mm
- 4499900 Size : Ø3 × L260 mm

### RET Control용 온도 센서

Cat. No.	Model	Type	재질	Diameter	EA Price ₩
H03-05-568	3736000	PT100.70	Glass Coated Stainless	Ø7 mm	540,000

- Length : 230 mm

### RCT basic & HS Digital용 온도 센서

Cat. No.	Model	Type	Description	EA Price ₩
H03-05-566	3516800	PT1000.60	Stainless Steel	300,000

- Size : Ø3 × L230 mm

### RCT Digital & HS Control용 온도 센서

Cat. No.	Model	Type	Description	EA Price ₩
H03-05-570	4443000	PT1000.80	Stainless Steel	370,000

- Size : Ø3 × L150 mm

### RET control-visc용 온도 센서, Stainless Steel

Cat. No.	Model	Type	Description	Size (mm)	EA Price ₩
H03-05-550	20000440	PT100.70	Single Probe	Ø3 × L230	300,000
H03-05-554	3367600	PT1000.50	Dual Probe		670,000
H03-05-542	2600300	PT100.51	Single Probe, Glass-coated	Ø8 × L230	690,000
H03-05-560	3377700	PT1000.51	Dual Probe, Glass-coated		770,000
H03-05-546	4499700	PT100.53	Single Probe, SafeCoat	Ø3 × L230	680,000

### RET control-visc Accessory

Cat. No.	Model	Description	내경(mm)	Max. Temp.	EA Price ₩
H03-05-691	20003416	HG 600 히팅자켓	91	Max. 300°C	380,000
H03-05-693	20003415	HG 1000 히팅자켓	106	Max. 300°C	510,000



ETS-D5



ETS-D6



H62.51



H70



H66.51



PT100.70



PT1000.60



PT1000.80



20003416



20003415



20007259

Cat. No.	Model	Description	작동시간	EA Price ₩
H03-05-699	20007259	Wireless Remote Control	15 h	970,000