

Digital Multi Hotplate Stirrer

디지털 멀티 가열 자력 교반기

Cat. No.	Model	Position	Plate Size	EA Price ₩
H03-06-300	HS20-M3	3	140×140 mm	980,000
H03-06-305	HS21-M3-Set	3		1,200,000
H03-06-310	HS20-M6	6 (3×2)		1,750,000

- 여러 개의 Plate를 사용 할 수 있어 공간 효율이 좋음
- 각각의 Plate는 독립적인 가열과 교반 제어가 가능함
- 내열성과 내화학성이 강하고 청소가 편리한 세라믹 코팅 상판적용
- 고사양 컨트롤러
 - 정확한 온도 제어와 RPM 제어가 가능한 Digital Feedback Control
 - RPM 설정이 편리한 조그 다이얼 스위치
 - 시안성이 좋고 각종 상태를 모니터링 할 수 있는 Back Light LCD Display
- 온도 제어와 과열방지 2가지 독립적인 센서 적용
- 내구성이 뛰어나고 소음이 적은 Brushless Motor 적용
- 외부센서를 통해 실제 시료 온도 제어 가능 (HS21-M3-Set)
- HS21-M3-Set 구성 : 본체, 외부센서, Support Kit (Support Rod, Holding Rod)
- 마그네틱바 이탈 방지 기능
- PL Insurance (제조물 배상 책임 보험)
- LK Labkorea



Cat. No.		H03-006-300	H03-06-305	H03-06-310
Model		HS20-M3	HS21-M3-Set	HS20-M6
외부온도센서 적용		불가능	가능	불가능
Plate	Position	3ea		6ea (3×2)
	Size	140×140 mm		
Temperature	Range	Max. 350°C		
	Resolution	Display 0.1°C / Control 0.1°C		
Stirring	Capacity (water)	20 L		
	Speed	80 ~ 1500 rpm		
	Resolution	5 rpm		
Material	Plate	Ceramic Coated Aluminum		
	Body	Powder Coated Aluminum		
Display		Digital LCD with Back-Light		
Timer		99 hr 59 min		
Safety		Over Temp. Protection		
Motor		Brushless Shade Motor		
Environmental Conditions		Temp. 5 ~ 40 °C, Humidity < 80 %		
Size	Plate	w500×d240×h110 mm	w500×d430×h130 mm	
Electric Supply	Power	1 Phase, 230Vac, 50/60Hz		
	Consumption	400 W×3 (6.5 A)	400 W×6 (10.9 A)	

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
H03-06-141	ET106	External Sensor, Ø3.2×L250 mm	62,500
H03-06-192	HSSR20	Support Rod	13,100
H03-06-193	HSRC20	Rod Clamp	8,700
H03-06-194	HSHR20	Holding Rod	18,100
H03-06-195	SKIT20	Support Kit	37,800

- Support Kit 구성품 : Support Rod, Boss Head Clamp, Holding Rod
- Optional Accessory는 HS21-M3-Set 모델에만 사용 가능

