

Digital Multi Hotplate Stirrer 디지털 멀티 가열 자력 교반기, 4 x 10L

Cat. No.	Model	온도	Speed	EA Price ₩
H03-238-665	MS-H340-S4	Max. 340 °C	200 - 1500 rpm	1,807,500


- PID Control 개별 Temperature & Speed Control
- ∅135mm Stainless Steel Ceramic Coating Plate
- Maintenance Free Brushless DC Motor
- DLAB

Model	MS-H340-S4
String Position	4
Max. Volume	10 L x 4개
Display	LCD
Temp. Range	Max. 340 °C (Overheating Protection 420 °C)
External Temp. Sensor	PT1000 (Accuracy : ±0.2°C)
Overheating Protection	420 °C
Plate Material	Ceramic Coated Stainless Steel
Speed Range	200 ~ 1500 rpm
IP Protection	IP21
Size, Plate/Overall	∅135mm, w698×d270×h128 mm
Power	200 ~ 240V, 50/60 Hz, 515W x 4

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
H03-238-577	18900294	PT1000 Sensor with Support Clamp Set	140,400

- Support Kit 구성품 : Support Rod, Boss Head Clamp, Holding Rod



Analog Multi Hotplate Stirrer 아날로그 멀티 가열 자력 교반기

Cat. No.	Model	Position	Plate Size (w×d)	EA Price ₩
S07-72-045	MS-1206P	6 (1열 6구)	120×120 mm	1,014,000
S07-72-050	MS-1526P	6 (2열 3구)	150×150 mm	1,037,000
S07-72-053	MS-1504P	4 (2열 2구)		754,000
S07-72-055	MS-1203P	3 (1열)	120×120 mm	546,000
S07-72-060	MS-1503P		150×150 mm	565,000

- Temp. Range : Max. 380°C
- 사용 중 마그네틱바가 이탈하거나 급회전 하지 않음
- 저속과 고속으로 다양하게 사용 가능
- 가열, 교반을 동시 또는 개별 컨트롤 가능
- 내열성, 내약품성, 내구성 등이 우수
- Mi-Sung



S07-72-050



S07-72-060

Cat. No.	S07-72-045	S07-72-050	S07-72-053	S07-72-055	S07-72-060
Model	HS12-06P	HS15-26P	HS15-24P	HS12-03P	HS15-03P
Stirring Position	6		4	3	
Temperature Range	Max. 380°C				
Speed Display / Range	Scale / 100 ~ 1500				
Max. Magnetic Bar	∅10×L100 mm				
Plate Size	w120×d120 mm	w150×d150 mm		w120×d120 mm	w150×d150 mm
Plate Material	Ceramic Coated				
Heat Output	480W×6ea	480W×4ea		480W×3ea	
Power	220V, 50/60 Hz				