

General - Purpose Hydrometer 표준 비중[부액]계, L300 mm (정밀형)

- 액체비중계 : 액체의 밀도를 측정하는 기구
- 표준비중(부액)계로 주로 사용 (Divi. 0.001 g/ml)
- LK Labkorea

Full Set with 19 pcs & 19본조 Case - set

Cat. No.	Model	Range (g/ml)	EA Price ₩
H04-13-100	HD099	0.700 ~ 1.850	626,500

Individual Hydrometer, 개당 비중계

Cat. No.	Model	Range (g/ml)	EA Price ₩
H04-13-112	HD001	0.700 ~ 0.760	28,200
H04-13-116	HD002	0.760 ~ 0.820	28,200
H04-13-120	HD003	0.820 ~ 0.880	28,200
H04-13-124	HD004	0.880 ~ 0.940	28,200
H04-13-128	HD005	0.940 ~ 1.000	28,200
H04-13-132	HD006	1.000 ~ 1.060	28,200
H04-13-136	HD007	1.060 ~ 1.120	28,200
H04-13-140	HD008	1.120 ~ 1.180	28,200
H04-13-144	HD009	1.180 ~ 1.240	28,200
H04-13-148	HD010	1.240 ~ 1.300	28,200
H04-13-152	HD011	1.300 ~ 1.360	28,200
H04-13-156	HD012	1.360 ~ 1.420	28,200
H04-13-160	HD013	1.420 ~ 1.480	28,200
H04-13-164	HD014	1.480 ~ 1.540	28,200
H04-13-168	HD015	1.540 ~ 1.600	28,200
H04-13-172	HD016	1.600 ~ 1.660	28,200
H04-13-176	HD017	1.660 ~ 1.720	28,200
H04-13-180	HD018	1.720 ~ 1.780	28,200
H04-13-184	HD019	1.780 ~ 1.850	28,200



General - Purpose Hydrometer 표준 비중[부액]계, L165 mm (일반형)

- 액체 비중계 : 액체의 밀도를 측정하는 기구
- 표준 비중(부액)계로 주로 사용 (Divi. 0.003 g/ml)
- LK Labkorea

Full Set with 19 pcs & 19본조 Case - set

Cat. No.	Model	Range (g/ml)	EA Price ₩
H04-13-200	HD199	0.700 ~ 1.850	109,600

Individual Hydrometer, 개당 비중계

Cat. No.	Model	Range (g/ml)	EA Price ₩
H04-13-212	HD101	0.700 ~ 0.760	4,700
H04-13-216	HD102	0.760 ~ 0.820	4,700
H04-13-220	HD103	0.820 ~ 0.880	4,700
H04-13-224	HD104	0.880 ~ 0.940	4,700
H04-13-228	HD105	0.940 ~ 1.000	4,700
H04-13-232	HD106	1.000 ~ 1.060	4,700
H04-13-236	HD107	1.060 ~ 1.120	4,700
H04-13-240	HD108	1.120 ~ 1.180	4,700
H04-13-244	HD109	1.180 ~ 1.240	4,700
H04-13-248	HD110	1.240 ~ 1.300	4,700
H04-13-252	HD111	1.300 ~ 1.360	4,700
H04-13-256	HD112	1.360 ~ 1.420	4,700
H04-13-260	HD113	1.420 ~ 1.480	4,700
H04-13-264	HD114	1.480 ~ 1.540	4,700
H04-13-268	HD115	1.540 ~ 1.600	4,700
H04-13-272	HD116	1.600 ~ 1.660	4,700
H04-13-276	HD117	1.660 ~ 1.720	4,700
H04-13-280	HD118	1.720 ~ 1.780	4,700
H04-13-284	HD119	1.780 ~ 1.850	4,700

