

Clock/ Thermometer/ Humidity 벽걸이형 시계, 온/습도 겸용

Cat. No.	Model	Humidity Range	EA Price ₩
H05-147-130	1071	0 ~ 100% RH	65,400

- Wall Mount Screw 기본 포함
- 국내 소급가능한 Traceable® 성적서 제공 (ILAC)
- Cole-Parmer / Traceable®

Model		1071
Humidity	Range	0 ~ 100% RH
	Resolution	2% RH
	Accuracy	±5% RH
Temperature	Range	-20 ~ 55°C
	Resolution	1°C
	Accuracy	±1.5°C
Clock Accuracy		0.01% (Traceable only)
Size		Ø254 mm
Power		AA battery



Clock / Thermometer / Humidity 벽걸이형 시계, 온/습도 겸용

Cat. No.	Model	Humidity Range	EA Price ₩
H05-147-173	1079	8 ~ 92% RH	156,000

- Wall Mount Screw 기본 포함
- 국내 소급가능한 Traceable® 성적서 제공 (ILAC)
- Cole-Parmer / Traceable®

Model		1071
Humidity	Range	8 ~ 92% RH
	Resolution	2% RH
	Accuracy	±5% RH
Temperature	Range	-37 ~ 47°C
	Resolution	2°C
	Accuracy	±1.5°C
Clock Accuracy		0.01% (Traceable only)
Size		w280 x d318 x h25.4 mm
Power		AA battery



Printing Humidity / Temperature Meter 프린트 내장형 온습도계

Cat. No.	Model	습도 범위	온도 범위	EA Price ₩
H05-147-050	4181	30 ~ 95% RH	-20 ~ 60°C / -200 ~ 1333°C	1,235,700

- Programmed interval 또는 버튼을 누르면 자동으로 습도, 온도 (Temperature1, Temperature2), 이슬점, 날짜, 시간이 프린트되어 제공
- 구성 : Main unit, Humidity/T1 temperature probe, T2 type K probe, Paper, Probe Holder/ Stand, ABS plastic case, AAA battery x 7ea
- 국내 소급가능한 Traceable® 성적서 제공 (ILAC)
- Cole-Parmer / Traceable®

Model		4181
Humidity	Range	30 ~ 95% RH
	Resolution	0.1% RH
	Accuracy	±3% RH (midrange) / ±5% (elsewhere)
Temperature1	Range	-20 ~ 60°C
	Resolution	0.1°C
Temperature2	Range	-200°C ~ 1333°C
	Resolution	0.1°C(-200 ~ 999.9°C), 1°C(-1000 ~ 1333°C)
	Accuracy	±0.2% of the reading plus 1.8°C
Response Time		less than 10 Seconds
Probe Size		Ø19×L159 mm, 1.22 m Cable(T1) / Ø1.5 mm, 1.22 m Cable(T2)
Size		w76×d39×h191 mm
Power		6×AAA battery, AC adapter(option)



Optional Accessory

Cat. No.	Model	Description	EA Price ₩
H05-147-550	8182	Printing Paper (pack of 2 rolls)	12,500
H05-147-560	8183	Adapter (for AC operation)	75,300