

## Table Top Shaking Incubator 테이블 탑 진탕 배양기

- 용도 : Flask, Petri Dish, Tube, Elisa Plate 등 샘플 혼합 및 배양
- 시간, 온도, 속도 표시 및 설정 가능
- Brushless DC Motor의 사용으로 저소음 및 뛰어난 안정성 보장
- 뚜껑을 열면 자동으로 작동이 멈추는 Running Safety System
- 자동 전원 복구, 과열 방지 장치 내장, Temperature Calibration 기능
- 다양한 Platform 장착 가능 (별도 구매)
- LK Labkorea

### Compact Table Top Shaking Incubator 230×230mm Platform

Cat. No.	Model	Temp. Range	Speed	EA Price ₩
S04-480-031	LS-0200	Ambient +5°C ~ 60°C	50 ~ 300 rpm	2,759,400

### Table Top Shaking Incubator 350×350mm Platform

Cat. No.	Model	Temp. Range	Speed	EA Price ₩
S04-480-035	LS-0300	Ambient +5°C ~ 60°C	50 ~ 250 rpm	3,608,400


Cat. No.	S04-480-031	S04-480-035
Model	LS-0200	LS-0300
Speed Range	50 ~ 300 rpm	50 ~ 250 rpm
Plate Size	230 x 230 mm	350 x 350 mm
Shaking Motion	Orbital	
Mixing Orbit	20 mm	
Controller	Control	Microprocessor Thermo Controller
	Display	LCD
	Timer	1 min ~ 99h 59min
Temperature	Range	Ambient +5°C ~ 60°C
	Accuracy	Max. ±0.5°C (at 37°C)
	Uniformity	Temperature Calibration Function
Heating Time	Max. 15 min (25 → 60°C)	Max. 20 min (25 → 60°C)
Internal Chamber Height	182mm	325mm
Size	w421 × d320 × h338 mm	w632 × d502 × h512 mm
무게	4 kg	10 kg
Electric Supply	Power	220 V, 50/60 Hz
	Fuse	250 V, 3A, Ø5×20
	Max Consumption	350 W      600 W

### Optional Platform for LS-0200

Cat. No.	Model	Description	EA Price ₩
S04-480-951	LS-1	for 5×250 ml Conical Flask ③	254,700
S04-480-954	LS-2	for 9×100 ml Conical Flask	254,700
S04-480-957	LS-3	for 10×50 ml Tube ②	212,300
S04-480-960	LS-4	for 18×15 ml Tube	212,300
S04-480-963	LS-5	Spring Wire Rack ①	254,700
S04-480-966	LS-6	for 4 Elisa Plate or Deep-Well Plate	254,700
S04-480-969	LS-7	for 5 Petri Dish	84,800

### Optional Platform for LS-0300

Cat. No.	Model	Description	EA Price ₩
S04-480-972	LS-T1	for 9×250 ml Conical Flask	424,500
S04-480-975	LS-T2	for 25×100 ml Conical Flask	424,500
S04-480-978	LS-T3	for 16 Petri Dish ④	148,700
S04-480-981	LS-T4	Spring Wire Rack	424,500



LS-0200



LS-0300

