

Grinding Ball

분쇄용 볼

Porcelain Ball / 자기 볼

Cat. No.	Model	Size	포장	EA Price ₩
M03-112-326	PB-20	Ø20 mm	1 kg	15,100
M03-112-332	PB-30	Ø30 mm		15,100
M03-112-338	PB-40	Ø40 mm		15,100
M03-112-344	PB-50	Ø50 mm		15,100
M03-112-350	PB-60	Ø60 mm		15,100
M03-112-356	PB-70	Ø70 mm		15,100

• LK Labkorea

Alumina Ball / 알루미나 볼

Cat. No.	Model	Size	포장	EA Price ₩
M03-110-352	AB93-01	Ø1 mm	1 kg	16,900
M03-110-360	AB93-02	Ø2 mm		16,900
M03-110-368	AB93-03	Ø3 mm		16,900
M03-110-384	AB93-05	Ø5 mm		16,900
M03-110-392	AB93-10	Ø10 mm		16,900
M03-110-408	AB93-20	Ø20 mm		16,900
M03-110-424	AB93-30	Ø30 mm		16,900

• Al₂O₃ (93 %), Material Density (3.6 g/cm³)

• LK Labkorea

Zirconia Ball / 지르코니아 볼 (Y₂O₃ZrO) 3mol% yttria (Y₂O₃) Stabilized Zirconia (ZrO₂)

Cat. No.	Model	Size	포장	EA Price ₩
M03-196-553	4039GM-S002	Ø0.2 mm	1 kg	300,700
M03-196-555	4039GM-S003	Ø0.3 mm		300,700
M03-196-557	4039GM-S005	Ø0.5 mm		156,600
M03-196-584	4039GM-S008	Ø0.8 mm		156,600
M03-196-593	4039GM-S010	Ø1 mm		144,600
M03-196-602	4039GM-S020	Ø2 mm		144,600
M03-196-611	4039GM-S030	Ø3 mm		144,600
M03-196-629	4039GM-S050	Ø5 mm		144,600
M03-196-634	4039GM-S065	Ø6.5 mm		144,600
M03-196-636	4039GM-S080	Ø8 mm		144,600
M03-196-638	4039GM-S100	Ø10 mm		144,600
M03-196-642	4039GM-S120	Ø12 mm		144,600
M03-196-647	4039GM-S150	Ø15 mm		144,600
M03-196-656	4039GM-S200	Ø20 mm		144,600
M03-196-666	4039GM-S250	Ø25 mm		144,600
M03-196-970	4039GM-S300	Ø30 mm		144,600

• 94.7wt% ZrO₂ + 5.3wt% Y₂O₃

• Material Density: 6.0 ~ 6.05 g/cm³

• Hardness: 11 GPa

• Inframat Advanced Material

Agate Ball / 아게이트 볼

Cat. No.	Model	Size	PK Price ₩
M03-81-502	604-K-002	Ø2 mm	10ea 8,800
M03-81-506	604-K-003	Ø3 mm	10ea 10,000
M03-81-510	604-K-005	Ø5 mm	10ea 14,300
M03-81-518	604-K-010	Ø10 mm	10ea 26,600
M03-81-526	604-K-015	Ø15 mm	10ea 85,700
M03-81-530	604-K-020	Ø20 mm	10ea 166,500
M03-81-534	604-K-025	Ø25 mm	10ea 283,500

• C.Giese



Porcelain Ball



Alumina Ball



Zirconia Ball



Agate Ball

B

Ball

Grinding Ball - Porcelain / Alumina / Zirconia / Agate