


Glass Sample Jar 유리 샘플 자 / 병, LukeGL®

Cat. No.	Model	용량	Size	Cap	EA Price ₩
J01-10-110	JK-010	100 ml	∅54×h98 mm	GL14	37,100
J01-10-117	JK-025	250 ml	∅70×h128 mm	GL18	47,200
J01-10-124	JK-050	500 ml	∅85×h148 mm	GL25	58,600

- Borosilicate Glass 3.3
- Teflon Seal이 내장된 Red PBT Cap (내화학적, 내열성)
- Strong Reagent 및 Solvent 등과 같이 사용 가능
- LukeGL®



Glass Specimen Jar 종자병, LukeGL®

Cat. No.	Model	용량	Size	Hole	EA Price ₩
J01-96-211	JD-010	100 ml	∅53×h98 mm	∅18 mm	2,800
J01-96-218	JD-025	250 ml	∅70×h128 mm	∅20 mm	3,800
J01-96-225	JD-050	500 ml	∅83×h148 mm	∅26 mm	5,900

- 국산 경질 유리
- Rubber Stopper 기본 포함
- LukeGL®



Glass Specimen Jar with Cover, Simax® 유리 스페시멘 자

Cat. No.	Model	Jar od	Cover od	Height	EA Price ₩
J01-985-111	416130808	85 mm	95 mm	85 mm	68,100
J01-985-115	416131010	100 mm	110 mm	100 mm	74,300
J01-985-119	416131212	120 mm	130 mm	120 mm	92,400
J01-985-123	416131515	150 mm	162 mm	150 mm	121,400
J01-985-127	416132020	200 mm	214 mm	200 mm	236,500
J01-985-131	416132525	250 mm	264 mm	250 mm	344,000
J01-985-135	416133030	300 mm	314 mm	300 mm	444,800

- Simax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- KAVALIERGLASS



Glass Jar, Rectangular, Sodalime Glass 직사각형 유리 자

Cat. No.	Model	Size	EA Price ₩
J01-83-101	30303.LK	w130×d170×h210 mm	83,000
J01-83-104	30304.LK	w140×d140×h220 mm	76,600
J01-83-107	30305.LK	w125×d225×h225 mm	134,200
J01-83-110	30306.LK	w150×d250×h250 mm	159,800

- Sodalime Glass
- LK Labkorea



더 많은 Glass Jar는 131 ~ 133 Page 참조

