

## Jar Tester 자 테스터

Cat. No.	Model	용량	교반 수	EA Price ₩
J03-03-130	MS-SF4	Max. 2 L	4개	845,000
J03-03-140	MS-SF6		6개	1,176,000

- 주로 물의 응집 실험을 위해 사용 되어짐
- 4 Place 용은 Portable Type으로 이동이 쉬움
- 하부와 후면에 고정된 조명장치는 강렬하지 않게 Sample로 조명을 확산시키므로 시료의 정확한 관찰이 용이함
- 내구성이 우수한 Brushless DC Motor를 사용 제품의 수명이 길고 안정적인임
- 속도조절은 Feedback 제어에 의해 Paddle의 회전속도를 정확하게 유지함
- 교반의 속도표시는 Digital Display되며, Timer 기능이 내장되어 있음 (99 hr 59 min)
- 설정 rpm과 실제 rpm이 선명한 LED화면에 각각 표시되어 동시에 확인이 가능
- Timer 설정에 의해서 정회전, 역회전, 자동 정·역회전을 선택할 수 있는 선택모드 기능
- Mi-sung

Cat. No.	J03-03-130	J03-03-140
Model	MS-SF4	MS-SF6
교반 위치 수	4 Place	6 Place
Stirring Capacity	Max. 2 L	
Motor	Brushless DC Motor 60 W	
Speed	20 ~ 300 rpm	
Speed Accu.	± 5 rpm at 20 ~ 49 rpm / ± 1 rpm at 50 ~ 300 rpm	
Controller	Feedback Control	
Rotation	Right / Left / Auto Reverse	
Timer	99 hr 59 min	
Size	w320×d320×h390 mm	w900×d250×h450 mm
Power	220V, 50/60 Hz	
Weight	9 kg	21.5 kg

- Auto Reverse : Reverse Time Setting or Time Operated



J03-03-130



J03-03-140

## Knife 나이프

Cat. No.	Model	길이	Blade	EA Price ₩
S02-139-221	1420-15	150 mm	45 mm	7,000
S02-139-231	1422-14	140 mm		6,000

- 재질 : 18/10 Stainless Steel, Antimagnetic, Chemical-Resistant
- LK Labkorea



1420-15



1422-14

## Laboratory Knife 실험실용 나이프

Cat. No.	Model	길이	Blade	EA Price ₩
S02-139-211	1416-18	180 mm	80 mm	8,000

- 재질 : 18/10 Stainless Steel, Antimagnetic, Chemical-Resistant
- LK Labkorea



## Scalpel Knife / Scalpel with Blade 메스 나이프

Cat. No.	Model	길이	Grip	EA Price ₩
S02-139-191	1412-15	150 mm	민무늬	4,200
S02-139-201	1414-15		빗살무늬	4,200

- 재질 : 18/10 Stainless Steel, Antimagnetic, Chemical-Resistant
- LK Labkorea



1412-15



1414-15