

## Portable Microscope 휴대용 현미경

Cat. No.	Model	Description	EA Price ₩
M02-29-402	MS-1	Portable Microscope	211,800

- 표면 측정을 위한 단안 현미경
- Eyepiece (WF10X / 15.5mm), 10X Achromatic Objective, 총배율 : 100X
- Penlight 부착형
- Optika Microscope

### Spare Penlight

Cat. No.	Model	Description	EA Price ₩
M02-29-999	M-899	Pen Illuminator	17,200
M02-29-912	15104	Cleaning Kit	24,800
M02-29-414	DC-001	Plastic dust cover, small, for B-60, MS/SFX & Ecovision Series 340(l)x400(h) mm	15,300



## Monoscope, Entry-level 입체 현미경, 교육용

Cat. No.	Model	Description	EA Price ₩
M02-29-412	MS-2	Monoscope	133,500

- Eyepiece (WF10X /16mm), 2X Achromatic Objective, 총배율 : 20X
- Working Distance : 57mm
- Head : 단안현미경, 45° 경사, 360° 회전가능
- Optika Microscope

### Optional Accessory

Cat. No.	Model	Description	EA Price ₩
M02-29-414	DC-001	Plastic dust cover, small, for B-60, MS/SFX & Ecovision Series 340(l) x 400(h) mm	15,300
M02-29-912	15104	Cleaning kit	24,800



## Compact, all-in-one Video Microscope 소형 비디오 현미경

Cat. No.	Model	Description	EA Price ₩
M02-29-432	IS-01	Video Microscope	3,149,600

- HDMI 카메라 및 Software 내장
- HDMI 케이블을 이용하여 모니터와 연결하여 사용 가능함, 현미경 접안렌즈를 통해 보는 것과 동일함
- 최대 50X 배율 까지 가능
- 전체 해상도에서 초당 60프레임 및 1080 HD 이미지 품질 제공
- Maximum Sample Height : 150 mm
- Objective Type : Achromatic
- Working Distance : 110 mm
- Standard Magnifications : 0.7x-5x
- Zoom Ratio : 7.14
- Zoom Click Stops\_ At 0.7x, 1x, 1.5x, 2x, 2.5x, 3x, 4x, 5x
- Full image Size : 1920x1080
- Frame Rate : 60 fps (1920x1080)
- Pixel Size : 2.9x2.9 μm
- Illuminance : 60000 lux
- Input Voltage : 110/240 Vac, 50/60 Hz
- Optika Microscope

