

## Monoscope Cordless, Entry-level 실험실 무선 생물현미경, 교육용

- 간단한 3단계 셋팅
  1. 선호하는 밝기를 선택
  2. 밝기 레벨 조절
  3. 렌즈를 변경하거나 조리개를 닫음 (현미경은 동일한 광도를 유지)
- 같은 렌즈로 오일과 물을 같이 사용 가능
  - 오일은 최상의 수치를 얻을 수 있는 최고의 미디어이며, 물은 적절한 결과와 편리함을 동시에 얻을 수 있다. 따라서 물은 교육용으로 사용하기를 권장
- LED Long Life Time (65,000시간)
- Head는 360° 회전이 가능하며, 30° 기울어져 있어 관찰이 용이함
- 눈과 눈 사이의 거리(Interpupillary)를 조절할 수 있어 누구나 쉽게 관찰할 수 있음 (Binocular)
- WF10X/18 mm 렌즈로 안경을 착용한 채로 관찰할 수 있음
- Stage : Fixed (고정), Mechanical (X,Y축 이동형)
- 구성품 : 본체 + 접안렌즈(Eyepiece) + 대물렌즈(Achromatic Objective)
- Optika Microscope

### Monocular Microscope

Cat. No.	Model	objective	Stage	EA Price ₩
M02-29-670	B-61	400x(4x,10x,40x)	Fixed stage	331,900
M02-29-672	B-62	400x(4x,10x,40x)	Mechanical stage	457,800
M02-29-674	B-63	600x(4x,10x,40x,60x)		541,800
M02-29-676	B-65	1000x(4x,10x,40x,100x)		541,800

### Binocular Microscope

Cat. No.	Model	objective	Stage	EA Price ₩
M02-29-678	B-66	400x(4x,10x,40x)	Mechanical stage	555,200
M02-29-680	B-67	600x(4x,10x,40x,60x)		641,000
M02-29-682	B-69	1000x(4x,10x,40x,100x)		641,000



Monocular



Binocular

Model	B-61	B-62	B-63	B-65	B-66	B-67	B-69
Eyepiece	Monocular(단안), WIDE FIELD 10X/18mm				Binocular(쌍안), WIDE FIELD 10X/18mm		
Head	Monocular, 360° rotating, 45° inclined				Binocular, 360° rotating, 30° inclined		
Stage	Fixed stage, 120x110 mm	Mechanical stage, 125x125mm, 62x24 mm(X-Y축이동)					
Condenser	N.A. 0.65 with iris diaphragm	Abbe N.A. 1.25, pre-centered, focusable, with iris diaphragm					
Illumination	1 W LED, with brightness control, rechargeable batteries. Color temperature: 6,300 K. Multi-plug 100-240Vac/5Vdc external power supply						

### Accessories

Cat. No.	Model	Description	EA Price ₩
M02-29-893	M-001	Huygens 5x eyepiece	45,800
M02-29-620	M-002.2	WF10x/18 eyepiece	36,300
M02-29-622	M-003.2	WF15x/12 eyepiece	36,300
M02-29-624	M-004.2	WF10x/18 eyepiece, with micrometric scale (10mm/100um)	64,800
M02-29-725	M-162	WF10x/10 eyepiece	57,200
M02-29-626	M-008.2	WF10x/18 eyepiece, with pointer	41,900
M02-29-628	M-131	Achromatic objective 4x/0.10	36,300
M02-29-630	M-132	Achromatic objective 10x/0.25	47,700
M02-29-632	M-133	Achromatic objective 20x/0.40	70,600
M02-29-634	M-134	Achromatic objective 40x/0.65	84,000
M02-29-636	M-135	Achromatic objective 60x/0.85	146,900
M02-29-638	M-136	Achromatic objective 100x/1.25 (oil)	150,700
M02-29-912	15104	Cleaning Kit	24,800



Cleaning Kit