

Digital Photo Video 디지털 사진 비디오

• Optika Microscope

USB Camera – Pro Series

Cat. No.	Model	Resolution	Frame Rate	EA Price ₩
M02-29-831	C-P3	2048×1536	50 frame / sec	1,201,800
M02-29-835	C-P6	3072×2048	30 frame / sec	1,442,300
M02-29-839	C-P8	3840×2160	32 frame / sec	1,499,500

Cat. No.	M02-29-831	M02-29-835	M02-29-839
Model	C-P3	C-P6	C-P8
Signal Output	USB 3.0		
Sensitivity	600mV(at 1/30 sec)	425mV(at 1/30 sec)	236mV(at 1/30 sec)
Sensor Type	Sony Exmor CMOS		
Pixel Size	2.5×2.5 um	2.4×2.4 um	1.62×1.62 um
Exposure Time	0.1m sec ~ 15 sec		
Power	PC USB		
C-Mount	Yes		

USB- Camera & Tablet PC – TB Serie

Cat. No.	Model	Resolution	Frame Rate	EA Price ₩
M02-29-841	TB-3W	2048×1536	12 frame / sec	1,728,400
M02-29-845	TB-5W	2592×1944	7 frame / sec	1,926,700

Cat. No.	M02-29-841	M02-29-845
Model	TB-3W	TB-5W
Signal Output	USB 2.0	
Sensitivity	1V/lux – second	0.53V/lux-second
Sensor Type	Aptina CMOS	
Pixel Size	3.2×3.2 um	2.2×2.2 um
Exposure Time	0.244 msec ~ 2 sec	0.294 msec ~ 2 sec
Power	PC USB	
C-Mount	Yes	
Tablet		
System	Window 10 64 bit	
CPU	Gemini-Lake, N4100	
CPU Speed	1.10 GHZ	
Graphic Card	Intel ® HD Graphics 600	
Memory	RAM 6 GB LPDDR3	
Display	LED 10.1" IPS Multi Touch Screen	
Resolution	1920×1200	

HDMI Camera – H Serie

Cat. No.	Model	Resolution	Frame Rate	EA Price ₩
M02-29-863	C-HB	1280×720	60 frame / sec	601,900
M02-29-876	C-HP4	3840×2160	30 frame / sec	1,864,900

Cat. No.	M02-29-863	M02-29-876
Model	C-HB	C-HP4
Signal Output	HDMI	HDMI, USB
Sensitivity	1300mV (at 1/30 sec)	505mV (at 1/30 sec)
Sensor Type	SONY STARVIS	SONY
Pixel Size	2.9×2.9 um	2.0×2.0 um
Exposure Time	0.01 m sec ~1 sec	0.04 m sec ~ 2 sec
Power	DC 12V 1A	
C-Mount	Yes	



다음 페이지에서 계속...