

Laboratory Stereomicroscopes

실험실 입체 현미경

- 기본 쌍안 (Binocular)이며, Photography/Video 기능을 장착할 수 있는 삼안(Trinocular) 모델 두종류임
- High Eyepoint 접안렌즈로 안경을 착용한 채로 관찰 할 수 있음
- Head는 360° 회전이 가능하며, 45° 기울어져 있어 관찰이 용이함
- Working Distance : 100 mm
- Optika Microscope

With Illuminator

Cat. No.	Model	Head	Objective	EA Price ₩
M02-29-125	SLX-1	Bionocular	Achromatic 2x...4x (with anti-fungus treatment)	611,800
M02-29-135	SLX-2	Bionocular	Achromatic 0.7x...4.5x (6.43 : 1 Ratio) with anti-fungus treatment	880,100
M02-29-145	SLX-3	Trinocular		974,900

- Cordless use, totally independent from mains/batteries connection
 - Rechargeable batteries, External Power Supply
- Longlife LED illumination

Without Illuminator

Cat. No.	Model	Head	Objective	EA Price ₩
M02-29-155	SLX-4	Bionocular	Achromatic 0.7x...4.5x (6.43 : 1 Ratio)	1,276,900
M02-29-165	SLX-5	Trinocular	with anti-fungus treatment	1,393,300

- Overhanging Stand : 420 mm high, 430 mm horizontal arm length, base 230x230 mm with following adjustments
- External illuminators available (Option)



SLX-1



SLX-4

Type	Binocular Type	Trinocular
Head	45° inclined, 360° rotating	Split Ratio 50/50, 45° inclined, 360° rotating
Dioptric adjustment	Both eyepieces	
Eyepiece	WF10x / 21 mm, high eye-point, Secured by screw and with rubber cups	
Working Distance	100 mm	
Focusing	Rack and pinion mechanism controlled by a pair of knobs placed on both side of the stand	

Optional Accessory

Eyepiece

Cat. No.	Model	Description	EA Price ₩
M02-29-201	ST-036	Eyepieces (pair), Flat	11,800
M02-29-203	ST-081	Eyepieces(pair) WF10x / 21mm, high eyepoint, w/ rubber cup	197,300
M02-29-205	ST-082	Eyepieces(pair) WF15x / 15mm, high eyepoint	197,300
M02-29-207	ST-083	Eyepieces(pair) WF20x / 10mm, high eyepoint	197,300
M02-29-209	ST-084	Eyepieces micrometric WF10x / 21mm, high eyepoint, w/ rubber cup	197,300

Camera Adapter

Cat. No.	Model	Description	EA Price ₩
M02-29-455	M-620	0.35x focusable C-mount adapter	351,000
M02-29-114	M-620.1	0.50x focusable C-mount adapter	360,600
M02-29-457	M-620.2	0.65x focusable C-mount adapter	400,600
M02-29-421	M-620.3	1.00x focusable C-mount adapter	144,900

Additional Lens

Cat. No.	Model	Description	EA Price ₩
M02-29-212	ST-085.1	Additional Lens 0.5x (w.d 165mm), SLX-1 제외	213,100
M02-29-214	ST-091	Additional Lens 0.75x (w.d 105mm), SLX-1 제외	193,400
M02-29-216	ST-086.1	Additional Lens 1.5x (w.d 45mm), w/ Compensating disc, SLX-1 제외	224,900
M02-29-218	ST-087	Additional Lens 2x (w.d 33mm), only for SLX-4 & 5	193,400

Additional Accessory

Cat. No.	Model	Description	EA Price ₩
M02-29-423	ST-088.1	Polarising Set	429,300
M02-29-426	ST-040.1	Darkfield Condenser	164,100
M02-29-427	ST-100.1	Hand Moving Stage	558,900
M02-29-556	CL-41	2-Arm LED Illuminator	665,100
M02-29-912	15104	Cleaning kit	24,800
M02-29-314	DC-002	Plastic Dust Cover for SLX 1,2,3	13,400
M02-29-316	DC-004	TNT Dust Cover for SLX 4,5	21,000