


## Slide Glass 슬라이드 글라스

Cat. No.	Model	Size	Type	PK Price ₩
M08-59-110	474743	w76×d26 mm	Plain ①	50ea 10,600
M08-59-120	474744		Twin Frosted ②	50ea 14,300

- 재질 : Optical Glass, Hydrolytic Class 3
- Complied with DIN ISO 8037-1
- Thickness : 1 mm
- Ground Edges
- Brand GmBH



## Slide Glass 슬라이드 글라스

Cat. No.	Model	Type	PK Price ₩	CS Price ₩
M08-660-111	0302-0002	Plain	50ea 4,100	20pk 76,500
M08-660-115	0302-2202	Double Frosted	50ea 4,700	20pk 86,900

- Size : w75×d25 mm, T1.1±0.1 mm
- Ground Edges, 45° Corners
- Pre-cleaned, Examined & Dried
- LK Labkorea



## Slide Glass 슬라이드 글라스

Cat. No.	Model	Size	Type	PK Price ₩
M08-625-110	LSG72P	w76×d26 mm	Plain ①	50ea 3,700
M08-625-120	LSG72F		Frosted ②	50ea 4,600

- Thickness : 1 mm
- LK Labkorea



## Color Frosted Slide Glass 컬러 슬라이드 글라스

Cat. No.	Model	Frosted Color	PK Price ₩	CS Price ₩
M08-660-123	0302-6101-02	Blue	50ea 6,100	20pk 112,900
M08-660-127	0302-6101-16	White	50ea 6,100	20pk 112,900
M08-660-131	0302-6101-17	Yellow	50ea 6,100	20pk 112,900

- Size : w75×d25 mm, T1.1±0.1 mm
- Ground Edges, 45° Corners
- Pre-cleaned, Examined & Dried
- LK Labkorea



## Large Slide Glass 대형 슬라이드 글라스

Cat. No.	Model	Size	PK Price ₩	CS Price ₩
M08-660-143	0351-0001	w76×d39 mm	100ea 16,200	5pk 77,200
M08-660-147	0350-0001	w76×d52 mm	100ea 18,300	5pk 87,000

- 두께 : 1.1±0.1 mm
- Grounded Edge, 90° Square Corners
- LK Labkorea

