

## Adhesion Slide Glass 코팅 슬라이드 글라스

Cat. No.	Model	Type	PK Price ₩	CS Price ₩
M08-660-119	0312-7141-16	White Frosted	50ea 15,900	20pk 296,400


- Permanent Especially-fomulated Adhesion으로 표면 코팅 : Hydrophobic Coating
- Tissue Section Loss나 Cell Loss의 위험성을 줄임
- Size : w75×d25 mm, T1.1±0.1 mm
- Ground Edges, 45° Corners
- Pre-cleaned, Examined & Dried
- LK Labkorea



## Cavity Slide Glass 홀 슬라이드 글라스

Cat. No.	Model	Size	Cavity수	PK Price ₩	CS Price ₩
M08-59-161	475505	w76×d26 mm	1	50ea 289,200	
M08-59-163	475535		2	50ea 434,300	
M08-59-165	475565		3	50ea 519,500	

- 재질 : Optical Glass, Hydrolytic Class 3
- Thickness : 1.2 ~ 1.5 mm
- Cavity : Ø15 ~ 18 mm, Depth : 0.6 ~ 0.8 mm
- Brand GmBH



## Cover Glass 커버 글라스

Cat. No.	Model	Type	Size (w×d)	PK Price ₩	CS Price ₩
M08-59-218	470045	정사각	18×18 mm	200ea 12,900	10pk 120,200
M08-59-226	470050		20×20 mm	200ea 17,500	10pk 163,500
M08-59-234	470055		22×22 mm	200ea 20,300	10pk 189,400
M08-59-238	470060	직사각	24×24 mm	200ea 24,600	10pk 229,800
M08-59-274	470816		24×40 mm	100ea 19,800	10pk 185,100
M08-59-286	470819		24×50 mm	100ea 26,200	10pk 244,200
M08-59-288	470820		24×60 mm	100ea 31,600	10pk 294,800

- Clear Borosilicate Glass, Hydrolytic Class 1
- Thickness : 0.13 ~ 0.17 mm, Flatness : ±3 um
- Retractive Index n<sub>e</sub> = 1.52±0.01
- Abbe Number V<sub>e</sub> = 56.5±0.5
- Dust & Grease Free, Distortion-free
- Brand GmBH



## Cover Glass 커버 글라스

Cat. No.	Model	Type	Size (mm)	PK Price ₩	CS Price ₩
M08-660-151	0330-0130	정사각	w18×d18	200ea 3,800	5pk 17,800
M08-660-155	0330-0630		w20×d20	200ea 3,900	5pk 20,100
M08-660-159	0330-1130		w22×d22	200ea 5,500	5pk 25,400
M08-660-163	0330-1630		w24×d24	200ea 5,900	5pk 27,600
M08-660-167	0330-2810	직사각	w24×d32	100ea 4,800	10pk 44,600
M08-660-171	0330-3110		w22×d40	100ea 5,200	10pk 47,800
M08-660-175	0330-3210		w24×d40	100ea 5,500	10pk 50,600
M08-660-179	0330-3610		w24×d50	100ea 6,500	10pk 60,100
M08-660-183	0330-4010		w24×d60	100ea 7,300	10pk 67,700

- 재질 : Borosilicate Glass 3.3
- Thickness : 0.13 ~ 0.16 mm • Uniform Flatness and Uniform Sizes
- Excellent Optical Properties, Stable Molecular Structure
- Hinged-lid Box in Individually Aluminum Vacuum Pack with Desiccant
- 포장 : Plastic Flip-flap Box (PK) / Aluminum Foil Tropical Wrapping (CS)
- LK Labkorea

