


Slide Glass Rack 슬라이드 글라스 랙

Cat. No.	Model	용량	EA Price ₩
M08-23-060	SGR-020	20매 (76×26 mm)	20,600


- 재질 : Stainless Steel 304
- Size : w81 × d150 × h22 mm
- 슬라이드 글라스 외 일반 Glass Plate의 스탠딩용
- Rayhope®



Slide Support 슬라이드 글라스 받침대

Cat. No.	Model	용량	EA Price ₩
R01-11-723	SGR-010	20매 (300×40 mm)	4,800


- 재질 : Stainless Steel Wire
- 슬라이드 글라스 받침대
- LK Labkorea



Slide Box 소형 슬라이드 박스, 5매용

Cat. No.	Model	거치 수	Size	EA Price ₩
M08-53-156	T482040	1매	w51 × d100 × h7 mm	1,300
M08-46-111	T482000	2매	w86 × d88 × h7 mm	1,100
M08-46-115	T482010	5매	w80 × d17 × h30 mm	1,200

- 재질 : PP, White Color
- 장식 개폐가 용이하게 제작됨
- Tarsons



Plastic Slide Mailer 사각 슬라이드 메일러

Cat. No.	Model	거치 수	EA Price ₩	PK Price ₩
M08-660-215	0500-1001	1매	200	200ea 36,500
M08-660-219	0500-2002	2매	270	100ea 26,000
M08-660-223	0500-1003	3매	400	100ea 36,500


- 재질 : PP, Natural Color
- LK Labkorea



Plastic Slide Storage Box 사각 슬라이드 박스

Cat. No.	Model	거치 수	EA Price ₩	PK Price ₩
M08-660-227	0500-1002	2매	300	100ea 19,400
M08-660-231	0500-1005	5매	200	100ea 15,200

- 재질 : PP, Natural Color
- LK Labkorea



Plastic Slide Tray 슬라이드 트레이

Cat. No.	Model	거치 수	EA Price ₩
M08-660-235	0500-5020-16	20매	4,400

- 재질 : PP, Natural Color
- LK Labkorea

