


Basket with Cover, PP 플라스틱 시험관망

Cat. No.	Model	Size (w×d×h)	EA Price ₩	PK Price ₩
B03-20-010	T180010	110×120×150 mm	7,800	6ea 44,700
B03-20-020	T180020	180×170×160 mm	8,300	6ea 47,300
B03-20-040	T180030	230×230×230 mm	12,700	6ea 72,800

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Test Tube 및 기타 Glassware의 보관 및 멸균용에 적합
- Cover는 탈부착이 가능
- Tarsons



PFA Basket & Tray PFA 바스켓 & 트레이

Cat. No.	Model	Type	내부 Size (w×d×h)	EA Price ₩
B03-64-111	F3200-002	Basket, Inner	273×207×140 mm	1,056,500
B03-64-115	F3200-001	Tray, Outer	275×210×143 mm	1,402,000

- 재질 : PFA, -196 ~ 250°C 내열
- 테프론 재질로 최상의 내화학성 가짐
- PFA Basket를 PFA Tray 안으로 겹쳐 놓을 수 있음
- 외부 Size : w365×d291×h152mm (Basket),
w374×d297×h147mm (Tray, 11L)
- Flon Industry



Stainless Steel Wire Basket, Square 스테인레스 정사각 시험관망

Cat. No.	Model	Size	EA Price ₩
B03-143-050	BSS150	w150×d150×h150 mm	25,000
B03-143-100	BSS180	w180×d180×h180 mm	26,300
B03-143-150	BSS200	w200×d200×h200 mm	29,500
B03-143-200	BSS250	w250×d250×h250 mm	44,700
B03-143-250	BSS300	w300×d300×h300 mm	61,800

- 재질 : Stainless Steel
- 소형 실험기구의 세척, 보관, 운반 및 멸균용에 적합
- LK Labkorea



Stainless Steel Wire Basket, Round 스테인레스 원형 시험관망

Cat. No.	Model	Size	EA Price ₩
B03-660-121	4052-0001	Ø150×h120 mm	24,700
B03-660-125	4052-0002	Ø210×h180 mm	37,600
B03-660-129	4052-0003	Ø240×h180 mm	39,400
B03-660-133	4052-0004	Ø300×h200 mm	54,900

- 재질 : Stainless Steel
- LK Labkorea



Stainless Steel Basket, Square 스테인레스 직사각 시험관망

Cat. No.	Model	Size	EA Price ₩
U02-56-817	LK-BA020	w120×d105×h65 mm	33,800
U02-56-823	LK-BA035	w215×d115×h65 mm	42,200
U02-56-829	LK-BA065	w274×d123×h96 mm	45,600
U02-56-835	LK-BA105	w262×d206×h93 mm	61,500
U02-56-841	LK-BA225	w452×d263×h93 mm	67,500

- 재질 : Stainless Steel
- LK Labkorea

