

Ball Mill Jar 볼밀용 포트

• LK Labkorea

Porcelain Mill Jar / 자재 포트 (자기 Material) ①

Cat. No.	Model	용량	Size	Type	EA Price ₩
M03-112-210	PMP-005	500 ml	Ø120×h150 mm	Stopper	84,300
M03-112-217	PMP-010	1 L	Ø150×h180 mm		113,000
M03-112-224	PMP-020	2 L	Ø180×h180 mm	Clamp	158,200
M03-112-231	PMP-030	3 L	Ø210×h240 mm		188,300
M03-112-238	PMP-050	5 L	Ø240×h270 mm		273,300
M03-112-245	PMP-100	10 L	Ø270×h300 mm		430,200

Alumina Mill Jar / 알루미나 포트 ②

Cat. No.	Model	용량	Size	EA Price ₩
M03-215-011	AMP-11	1.1 L	Ø140×h160 mm	1,018,100
M03-215-020	AMP-20	2.0 L	Ø165×h200 mm	1,245,700
M03-215-036	AMP-36	3.6 L	Ø180×h261 mm	1,500,400
M03-215-050	AMP-50	5.0 L	Ø210×h268 mm	2,229,200
M03-215-073	AMP-73	7.3 L	Ø240×h289 mm	3,043,700

Alumina Mill Jar / 알루미나 포트 ③

Cat. No.	Model	용량	Size	EA Price ₩
M03-110-052	AMP-01	500 ml	Ø100×h110 mm	833,300
M03-110-060	AMP-05	1 L	Ø123×h124 mm	952,300
M03-110-068	AMP-10	2 L	Ø154×h166 mm	1,522,400
M03-110-076	AMP-20	3 L	Ø170×h180 mm	2,335,300
M03-110-084	AMP-30	5 L	Ø200×h220 mm	3,414,500

Zirconia Mill Jar / 지르코니아 포트 ④

Cat. No.	Model	용량	Size	EA Price ₩
M03-196-161	ZMP-05	500 ml	Ø95×h110 mm	2,443,400
M03-196-169	ZMP-10	1 L	Ø122×h130 mm	3,649,400

Agate Mill Pot with Lid ⑤

Cat. No.	Model	용량	Size	EA Price ₩
M03-81-351	602-Z-010	10 ml	Ø100×h110 mm	870,500
M03-81-354	602-Z-050	50 ml	Ø123×h124 mm	732,800
M03-81-359	602-Z-080	80 ml	Ø122×h130 mm	1,008,600
M03-81-363	602-Z-250	250 ml	Ø154×h166 mm	1,901,700
M03-81-367	602-Z-500	500 ml	Ø170×h180 mm	4,252,200

- 고품질의 아게이트 분쇄용 포트
- PTFE Gasket 포함
- C.Giese

Stainless Steel Mill Jar / 스테인레스 포트 ⑥

Cat. No.	Model	용량	Size	EA Price ₩
M03-215-101	1-3951-01	430 ml	Ø89.1×h86 mm	575,100
M03-215-104	1-3951-02	1 L	Ø114.3×h116 mm	696,000
M03-215-107	1-3951-03	3 L	Ø163×h200 mm	877,500

- 재질 : Stainless Steel 304
- Silicone Rubber 포함
- Universal



Grinding Ball / 분쇄용 볼

- Porcelain, Alumina, Zirconia, Agate, Quartz, Glass Ball은
66 ~ 68 Page 참조

